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台塑白蛋白試劑（ALB）- Bromcresol green method

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**MeDiPro ALBUMIN TEST (ALB) - Bromcresol green method****INTENDED USE**

For the quantitative determination of albumin in serum or plasma.


CLINICAL SIGNIFICANCE

Albumin is the most abundant individual protein in serum. The liver is the site of albumin synthesis, hence liver disease is one of the common causes of decreased albumin level. Other conditions such as congestive heart failure, malnutrition, starvation, nephritis, gloemlonephritis, leukemia and severe burns also results in the decrease of albumin level.

Albumin levels may be elevated in dehydration, shock, hemoconcent ration, or upon administration of large quantities of concentrated albumin intravenously.

PRINCIPLE

Bromcresol green (BCG), an ionic dyes, binds tightly to albumin when serum is added. The albumin- BCG complex absorbs light at 600 nm more intensively than BCG alone. The reaction formula is as follows:

**REAGENT**

- Package: please see the reagent box label shown.
- Components:

Component	Conc.
Succinate buffer(pH4.2±0.2)	75 mmol/L
Bromcresol green	0.15 mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C.**, reagent should be protected from light and contamination should be avoided. **Do not freeze the reagent!**

SPECIMEN COLLECTION AND PREPARATION

Fresh serum, heparin or EDTA-treated plasma, and avoid hemolysis.

PROCEDURES

- Main wavelength : 600 nm
Sub. wavelength : 700nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains single reagent, ready to use.

	Blank	Control	Sample
Sample (ml)	---	---	0.01
Control (ml)	---	0.01	---
Dist. Water (ml)	0.01	---	---
Reagent (ml)	1.0	1.0	1.0

Mix, incubate at 37°C for 5 min, and read the absorbance against reagent blank. $A = A_{600} - A_{700}$

CALCULATION

With standard or calibrator

$$\text{Albumin (g/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (g/dL)}$$

REFERENCE RANGE

Age	Conc. (g/dL)
Adult	3.5-5.0
>60yr.	3.4-4.8
Pregnancy	2.5-4.0

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition.
- The test is developed to determine albumin concentrations up to 5.5 g/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)
- According to the technical suggestion, the volume of reagent and specimen could be

**MeDiPro ALBUMIN TEST (ALB) - Bromcresol green method**

adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[g/dL]	SD [g/dL]	CV[%]
Sample1	2.77	0.01	0.52
Sample2	2.80	0.01	0.51
Sample3	3.99	0.03	0.64

- Linearity

This kit has a good linearity up to 7.0g/dL.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 300mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 20mg/dL
Intrafat	No interference was observed by intrafat up to 2.0%

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

Doumas BT, Watson WA, Biggs HG. Albumin standards and the measurement of serum albumin with bromcresol green. Clin. Chim. Acta. 1971. 31:87-96.

PARAMETER SETUP**Hitachi 7170/917 Applications**

TEST	[ALB]
ASSAY CODE	[1POINT]:[10]-[0]
SAMPLE VOLUME	[2]
R1 VOLUME	[200]
R2 VOLUME	[0]
WAVELENGTH(nm)	[700][600]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[ALB]
ASSAY CODE	[1POINT]:[15]-[0]
SAMPLE VOLUME	[3]
R1 VOLUME	[300]
R2 VOLUME	[0]
WAVELENGTH(nm)	[700][600]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0005M	MeDiPro ALBUMIN TEST	R1 6x20ml
BC-0005A	MeDiPro ALBUMIN TEST	R1 4x60ml
BC-0005B	MeDiPro ALBUMIN TEST	R1 4x100ml
BC-0005C	MeDiPro ALBUMIN TEST	R1 4x300ml
BC-0005D	MeDiPro ALBUMIN TEST	R1 4x500ml

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan