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
**MeDiPro ALKALINE PHOSPHATASE ( ALP ) - IFCC method****INTENDED USE**

For the quantitative determination of alkaline phosphatase activity in serum.

**CLINICAL SIGNIFICANCE**

Alkaline phosphatases(ALP) in serum is derived primarily from the liver, intestine and with little amount from bone. Increase of ALP activities are, therefore, indicators of hepatobiliary or bone disorder. In hepatobiliary disease, biliary cirrhosis produces exceptionally high levels of serum ALP. Infiltrative disease and cholangiolitic hepatitis are accompanied by significant elevations of ALP. Widely elevated ALP level may also occur in infectious mononucleosis.

Among the bone diseases, the highest level of serum ALP activity happens in Paget's disease and bone cancer. Moderate rises are observed in osteomalacia, rickets, Faconi syndrome, primary hyperparathyroidism and secondary hyperparathyroidism. An increase of alkaline phosphatase up to 2 to 3 times of normal is observed in women in the third trimester of pregnancy, this enzyme is placental original. Moderate elevations of ALP may be observed in several disorders that do not involve the liver or bone. Among these are Hodgkin's disease, congestive heart failure, ulcerative colitis, regional enteritis, and intraabdominal bacterial infections.

**PRINCIPLE****REAGENT**

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R1	AMP (buffer) pH10.2	413 mmol/L
	MgCl <sub>2</sub>	1.25 mmol/L
R2	p-nitrophenyl phosphate	32.5 mmol/L

**STORE TEMPERATURE**

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C**., reagent should be protected from light and contamination should be avoided.

**Do not freeze the reagent!**

**SPECIMEN COLLECTION AND PREPARATION**

Freshly drawn serum is the specimen of choices. Avoid stasis or hemolysis. Heparinized plasma may also be used. Fluoride and anticoagulants such as citrate, oxalate or EDTA should be avoided in the collection of specimens due to their ability to be chelated with Mg<sup>2+</sup>. Alkaline phosphatase in serum is stable for at least 7 days at 4°C and longer when frozen.

**PROCEDURES**

- Main wavelength : 405 nm  
Sub. wavelength : 660nm  
Reaction Temperature : 37°C  
Optical path length : 1.0 cm

- This kit contains two reagents, ready to use.

	Volume (ml)
Sample	0.02
R <sub>1</sub>	0.8
Mix, 37°C incubate 5 min	
R <sub>2</sub>	0.2

Mix, incubate at 37°C for 1 min, and read the initial absorbance A<sub>1</sub> against reagent blank, then read end absorbance A<sub>2</sub> in every 1 min. A = A<sub>405</sub> - A<sub>660</sub>.

**CALCULATION**

$$\text{ALP (U/L)} = \frac{(A_2 - A_1) / \text{min} \times V_t \times 1000}{L_p \times \epsilon \times V_s}$$
$$= (A_2 - A_1) / \text{min} \times 2742$$

V<sub>t</sub>: Reaction total volume 1.02 ml, V<sub>s</sub>: sample volume 0.02 ml  
ε: p-nitrophenyl phosphate molar absorptivity 18.6,  
1000: transfer U/ml to U/L, L<sub>p</sub>: Optical path length (cm)

**REFERENCE RANGE**

Age	Conc. (U/L)
Child	30-130
Adult	20-70

**WARNINGS AND PRECAUTIONS**

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Red cells contains lots of alkaline phosphatase than serum, if specimen hemolyzed, it will interfere with the result.
- The test is developed to determine alkaline phosphatases concentrations up to 800U/L. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- Reagent R2 shall be clear light yellow solution. Please don't use when it becomes muddy or blank OD over 0.7.

**MeDiPro ALKALINE PHOSPHATASE ( ALP ) - IFCC method**

- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.

- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.

- Waste management please refers to the local legal requirements.

- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)

- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


- Validity please see the reagent box label shown.

**REAGENT CHARACTERS**

- Precision (Within run)

N=15	Mean[U/L]	SD [U/L]	CV[%]
Sample1	85.80	0.63	0.74
Sample2	378.8	3.61	0.95
Sample3	382.1	3.11	0.81

- Linearity

This kit has a good linearity up to 900U/L.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 150mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 2.0%

- Stability

Expire day	1 year
Open vial stability	14 day

**REFERENCE**

- Posen S, Doherty E. The measurement of serum alkaline phosphatase in clinical medicine. Adv. Clin. Chem. 1981. 157: 1441.
- Fleisher GA, Eickelberg ES, Elveback LR, Alkaline phosphatase activity in the plasma of children and adolescents. Clin. Chem. 1977. 23: 469.

**PARAMETER SETUP****Hitachi 7170/917 Applications**

TEST	[ALP]
ASSAY CODE	[Rate A]: [19]-[34]
SAMPLE VOLUME	[4]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH (nm)	[660][405]
CALIB. METHOD	[Linear]

**Hitachi 7150/717 Applications**

TEST	[ALP]
ASSAY CODE	[Rate A]: [30]-[50]
SAMPLE VOLUME	[6]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH (nm)	[660][405]
CALIB. METHOD	[Linear]

**ORDERING INFORMATION**

Cat. No.	Product	Package
BC-0006M	MeDiPro ALKALINE PHOSPHATASE TEST	R1 6×20ml R2 3×10ml
BC-0006A	MeDiPro ALKALINE PHOSPHATASE TEST	R1 4×60ml R2 2×30ml
BC-0006B	MeDiPro ALKALINE PHOSPHATASE TEST	R1 4×100ml R2 2×50ml
BC-0006C	MeDiPro ALKALINE PHOSPHATASE TEST R1	R1 2×300ml
BC-0006D	MeDiPro ALKALINE PHOSPHATASE TEST R1	R1 2×500ml
BC-0006G	MeDiPro ALKALINE PHOSPHATASE TEST R2	R2 2×200ml

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan