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**MeDiPro ASPARTATE AMINOTRANSFERASE (AST) - UV-IFCC method****INTENDED USE**

For the quantitative determination of aspartate aminotransferase activity in serum.

CLINICAL SIGNIFICANCE

Serum aspartate aminotransferase (AST) catalyzes the transfer of the amino group from aspartic acid to α-ketoglutaric acid. This enzyme is found in practically every tissue of the body, including red blood cells. The concentration is particularly high in the liver, heart and skeletal muscles. Acute destruction of tissue results in the release of AST into the blood stream. In myocardial infarction, there is a significant increase in serum AST activity in 6 to 8 hours with a peak value reached after 48 to 60 hours. However, serum ALT activity remains within normal limits or only marginally increased. In hepatitis and other forms of liver disease associated with hepatic necrosis, both AST and ALT are elevated. Elevated levels of serum AST activity are also observed in infectious mononucleosis, muscular dystrophy, dermatomyositis, and in other forms of muscle and liver injury. The method presented here is an UV-Kinetic method based on the rate of NADH oxidation in a coupled malic dehydrogenase reaction.

PRINCIPLE

L-aspartate + α-ketoglutaric acid $\xrightarrow{\text{AST}}$ oxalacetate + L-glutamate

Oxalacetate + NADH + H⁺ $\xrightarrow{\text{MDH}}$ L-malate + NAD⁺

REAGENT

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R ₁	Tris Buffer	80 mmol/L
	MDH	> 600 U/L
	L-aspartate	240 mmol/L
	α-ketoglutaric acid	12 mmol/L
R ₂	Tris Buffer	80 mmol/L
	NADH	0.18 mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at 2 – 8 °C., reagent should be protected from light and contamination should be avoided.

Do not freeze the reagent!

SPECIMEN COLLECTION AND PREPARATION

Either serum or plasma may be used. The use of oxalate, citrate, EDTA, or heparin has shown no effect on AST values. But serum is recommended to be the choice. Hemolysis must be avoided because AST activity in red cells is 10~40 times higher than that of plasma. The serum or plasma should be removed from the clot or cells without delay. AST is reported to be stable for 3 to 4 days at room temperature; 2 weeks when stored refrigerated at 2~8°C and longer when frozen.

PROCEDURES

- Main wavelength : 340 nm
Sub. wavelength : 405nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains two reagents, ready to use.

	Volume (ml)
Sample	0.05
R ₁	0.8
Mix, 37°C incubate 5min	
R ₂	0.2

Mix, incubate at 37°C for 1 min, and read the initial absorbance A₁ against reagent blank, then read end absorbance A₂ in every 1 min. A = A₃₄₀-A₄₀₅.

CALCULATION

$$\text{AST (U/L)} = \frac{(A_2 - A_1) / \text{min} \times V_t \times 1000}{L_p \times \epsilon \times V_s}$$
$$= (A_2 - A_1) / \text{min} \times 3376$$

V_t: Reaction total volume 1.05 ml, V_s: sample volume 0.05 ml

ε: NADH molar absorptivity 6.22,

1000: transfer U/ml to U/L, L_p: Optical path length (cm)

REFERENCE RANGE

<46U/L

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition.
- Do not use when reagent blank OD less than 1.0.
- The test is developed to determine aspartate aminotransferase concentrations up to 800U/L. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.

**MeDiPro ASPARTATE AMINOTRANSFERASE (AST) - UV-IFCC method**

- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.

- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.

- Waste management please refers to the local legal requirements.

- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)

- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[U/L]	SD [U/L]	CV[%]
Sample1	35.13	0.52	1.47
Sample2	181.40	0.91	0.50
Sample3	165.47	0.74	0.45

- Linearity

This kit has a good linearity up to 1000U/L.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 150mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 24mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 2.0%

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C, Part 5. Clin. Chem. Lab. Med. 2002, 40: 631.
- New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clinica Chimica. Acta. 2003, 327: 69.

PARAMETER SETUP**Hitachi 7170/917 Applications**

TEST	[AST]
ASSAY CODE	[Rate A]: [19]-[34]
SAMPLE VOLUME	[10]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH(nm)	[405][340]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[AST]
ASSAY CODE	[Rate A]: [30]-[50]
SAMPLE VOLUME	[15]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH(nm)	[405][340]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0008M	MeDiPro ASPARTATE	R1 6×20ml
	AMINOTRANSFERASE TEST	R2 3×10ml
BC-0008A	MeDiPro ASPARTATE	R1 4×60ml
	AMINOTRANSFERASE TEST	R2 2×30ml
BC-0008B	MeDiPro ASPARTATE	R1 4×100ml
	AMINOTRANSFERASE TEST	R2 2×50ml
BC-0008C	MeDiPro ASPARTATE	R1 4×300ml
	AMINOTRANSFERASE TEST R1	
BC-0008D	MeDiPro ASPARTATE	R1 4×500ml
	AMINOTRANSFERASE TEST R1	
BC-0008G	MeDiPro ASPARTATE	R2 4×200ml
	AMINOTRANSFERASE TEST R2	

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan