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台塑鈣試劑（CA）- Arsenazo III method

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MeDiPro CALCIUM TEST (CA) - Arsenazo III method

INTENDED USE

For the quantitative determination of calcium in serum or plasma.


CLINICAL SIGNIFICANCE

The metabolism of calcium and that of phosphorus are very closely related. More than 99% of the calcium in the body are present in bone as calcium phosphate, the remains of the calcium, although in very small amount, has various and significant functions in the body. For example, calcium ions will decrease neuromuscular excitability, participate in blood coagulation, and activate some enzymes, such as succinate dehydrogenase and adenosine triphosphatase.

Hypercalcemia occurs in hyperparathyroidism, hyper-vitaminosis, sarcoidosis, multiple myeloma and certain cancers of bone. Hypocalcemia is observed in hypothyroidism, rickets, nephrosis, nephritis, steatorrhea and pancreatitis. Serum calcium level, is slightly dependent upon protein concentration, especially albumin. Hypoalbuminemia is usually accompanied by the decreased serum calcium level.

PRINCIPLE

Calcium, reacts with Arsenazo to form a bluish-purple colored chromophore which with an absorbance at 650 nm. Magnesium and iron are excluded from the reaction by complexing with 8-hydroxyquinoline. The reaction equation is shown as following:

REAGENT

- Package: please see the reagent box label shown.
- Components:

Component	Conc.
Imidazole	28 mmol/L
8-HQ	2 mmol/L
Arsenazo-III	0.2 mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C.**, contamination should be avoided.

Do not freeze the reagent!

SPECIMEN COLLECTION AND PREPARATION

Freshly drawn serum is the specimen of choice, plasma derived from specimen collection tubes containing EDTA, citrate, or oxalate must not be used in this test. Serum should be removed from the clot without delay to avoid loss of serum calcium due to increased erythrocyte permeability to calcium. Serum calcium separated from clot is stable at least 24 hours at room temperature, 1 week under refrigeration and up to 1 year in the freezer.

PROCEDURES

- Main wavelength : 660 nm
Sub. wavelength : 700nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains single reagent, ready to use.

	Blank	Control	Specimen
ddH ₂ O (ml)	0.01	---	---
Control (ml)	---	0.01	---
Specimen (ml)	---	---	0.01
R1 (ml)	1.0	1.0	1.0

Mix, incubate at 37°C for 5 min, and read the absorbance against reagent blank. $A = A_{660} - A_{700}$

CALCULATION

With standard or calibrator

$$\text{Calcium (mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

REFERENCE RANGE

8.7-11.0 mg/dL (2.18-2.74 mmol/L)

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition.
- The instrument, material and distilled water shall avoid calcium contamination.
- The test is developed to determine calcium concentrations up to 14mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal requirements.

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

MeDiPro CALCIUM TEST (CA) - Arsenazo III method

- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)

- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	8.27	0.10	1.20
Sample2	11.64	0.13	1.13
Sample3	11.87	0.11	0.96

- Linearity

This kit has a good linearity up to 15mg/dL.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 500mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 32mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 1.8%

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

- Janssen J.W. and Helbing A.R. Arsenazo III: An improvement of the routine calcium determination in serum. Eur. J. Clin Chem. Clin. Biochem. 1991, 29:197-201.
- Leary N.O., Pembroke A., Duggan P. F. Single stable reagent (Arsenazo III) for optically robust measurement of calcium in serum and plasma. Clin. Chem. 1992, 38: 904-908

PARAMETER SETUP

Hitachi 7170/917 Applications

TEST	[CA]
ASSAY CODE	[1 POINT]:[6]-[0]
SAMPLE VOLUME	[2]
R1 VOLUME	[200]
R2 VOLUME	[0]
WAVELENGTH(nm)	[700][660]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[CA]
ASSAY CODE	[1POINT]:[10]-[0]
SAMPLE VOLUME	[3]
R1 VOLUME	[300]
R2 VOLUME	[0]
WAVELENGTH(nm)	[700][660]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0013M	MeDiPro CALCIUM TEST	R1 6x20ml
BC-0013A	MeDiPro CALCIUM TEST	R1 4x60ml
BC-0013B	MeDiPro CALCIUM TEST	R1 4x100ml
BC-0013C	MeDiPro CALCIUM TEST	R1 3x250ml
BC-0013D	MeDiPro CALCIUM TEST	R1 2x500ml

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan