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台塑膽固醇試劑（CHOL） - Enzymatic method

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MeDiPro CHOLESTEROL TEST (CHOL) - Enzymatic method


INTENDED USE

For the use of quantitative determination of cholesterol in serum or plasma.

CLINICAL SIGNIFICANCE

Total cholesterol in the blood is elevated in idiopathic hypercholesterolemia, primary and secondary hyperlipidemias, diabetes mellitus, nephrotic syndrome, hypothyroidism and biliary obstruction. Pregnancy may also accompanied by a moderate increase of cholesterol level. Decrease of cholesterol is seen in patients with severe hepatitis and occasionally in severe anemia or infection.

PRINCIPLE

REAGENT

- Package: please see the reagent box label shown.
- Components:

Component	Conc.
Phenol	10 mmol/L
Sodium cholate	50 mmol/L
4-AA	1.4 mmol/L
Cholesterol esterase	200 U/L
Cholesterol oxidase	180 U/L
peroxidase	3000 U/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at 2 – 8 °C., reagent should be protected from light and contamination should be avoided.

Do not freeze the reagent!

SPECIMEN COLLECTION AND PREPARATION

Serum or heparinized plasma is recommended. The use of oxalate, citrate, fluoride and EDTA might result in slightly lower cholesterol values.

PROCEDURES

- Main wavelength : 505 nm
Sub. wavelength : 660nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains single reagent, ready to use.

	Control	Specimen	Blank
R ₁ (ml)	1.0	1.0	1.0
Control (ml)	0.01	----	----
Specimen (ml)	----	0.01	----
ddH ₂ O (ml)	----	----	0.01

Mix, incubate at 37°C for 10 min, and read the absorbance against reagent blank. $A = A_{505} - A_{660}$

CALCULATION

With standard or calibrator

$$\text{Cholesterol (mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

REFERENCE RANGE

150~250 mg/dL

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Specimen should avoid ascorbic acid which might interfere with the test result.
- The test is developed to determine cholesterol concentrations up to 400mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)
- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.

MeDiPro CHOLESTEROL TEST (CHOL) - Enzymatic method


11. Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	243	2.50	1.03
Sample2	128	0.96	0.75
Sample3	131	1.08	0.83

- Linearity

This kit has a good linearity up to 500mg/dL.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 200mg/dL
Ascorbic acid	Not suitable when sample contain reducing agent
Bilirubin (free form)	No interference was observed by bilirubin up to 12mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 16mg/dL
Intrafat	No interference was observed by intrafat up to 1.0%

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

- Allain, cc, Poon, LS, Chan, LS, et.al. Enzymatic determination of total serum cholesterol. Clin. Chem. 1974, 20:470.
- Bishop ML, Duben-Engelkirk JL, Fody EP. Clinical Chemistry: principles, procedures, correlations. Fourth edition, 2000.

PARAMETER SETUP

Hitachi 7170/917 Applications

TEST	[CHOL]
ASSAY CODE	[1 POINT]: [34]-[0]
SAMPLE VOLUME	[2]
R1 VOLUME	[200]
R2 VOLUME	[0]
WAVELENGTH (nm)	[660][505]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[CHOL]
ASSAY CODE	[1 POINT]: [50]-[0]
SAMPLE VOLUME	[2]
R1 VOLUME	[200]
R2 VOLUME	[0]
WAVELENGTH(nm)	[660][505]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0014M	MeDiPro CHOLESTEROL TEST	R1 6×20ml
BC-0014A	MeDiPro CHOLESTEROL TEST	R1 4×60ml
BC-0014B	MeDiPro CHOLESTEROL TEST	R1 4×100ml
BC-0014C	MeDiPro CHOLESTEROL TEST	R1 4×300ml
BC-0014D	MeDiPro CHOLESTEROL TEST	R1 4×500ml

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan