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台塑肌酸激酶試劑（CK）- UV kinetic, NAC method

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MeDiPro CREATINE KINASE TEST (CK) - UV kinetic, NAC method

INTENDED USE

For the quantitative determination of creatine kinase activity in serum or plasma.

CLINICAL SIGNIFICANCE

Creatine kinase (CK) catalyzes the reversible reaction of creatine and ATP to form creatine phosphate and ADP which plays a very important role in the function of energy storage in the human tissue. It presents mostly in skeletal muscle, heart and brain. Damage of muscle and heart tissues will release CK to the blood stream and result in the increase of CK activity. Determination of CK activity in serum or plasma is a good diagnostic marker of muscular dystrophy, myocardial infarction, renal damage and dysfunction and other skeletal muscle diseases.

PRINCIPLE

Creatine phosphate + ADP $\xrightarrow{\text{CK}}$ creatine + ATPATP + D-glucose $\xrightarrow{\text{Hexokinase}}$ glucose-6-phosphate + ADPglucose-6-phosphate + NAD⁺ $\xrightarrow{\text{G6PDH}}$ 6-phosphogluconate + NADH + H⁺

The rate of NADH formation is directly proportional to the CK activity.

REAGENT

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R ₁	Imidazole pH7.0±0.3	10 mmol/L
	Mg ²⁺	10 mmol/L
	AMP	5.0 mmol/L
	Hexokinase	3000 U/L
	G6PDH	2500 U/L
R ₂	Creatine phosphate	30 mmol/L
	D-glucose	20 mmol/L
	ADP	5 mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at 2 – 8 °C., contamination should be avoided.

Do not freeze the reagent!

SPECIMEN COLLECTION AND PREPARATION

Serum and heparin-treated plasma are the choices for the assay. Avoid from hemolysis since glucose-6-phosphate dehydrogenase, adenylate kinase, ATP liberated from red cells will interfere with the result. CK is reportedly stable for 4 hrs at room temperature, 8~12 hours at 4°C., and 2~3 days when frozen. Analyze as soon as possible after sampling is recommended.

PROCEDURES

- Main wavelength : 340 nm
Sub. wavelength : 405nm
Reaction Temperature : 37°C

Optical path length : 1.0 cm

- This kit contains two reagents, ready to use.

	Volume (ml)
sample	0.02
R ₁	0.8
Mix, 37°C incubate 3min	
R ₂	0.2

Mix, incubate at 37°C for 3 min, and read the initial absorbance A₁ against reagent blank, then read end absorbance A₂ in every 1 min. A=A₃₄₀-A₄₀₅.

CALCULATION

$$\text{CK (U/L)} = \frac{(A_2 - A_1) / \text{min} \times V_t \times 1000}{L_p \times \epsilon \times V_s}$$

$$= (A_2 - A_1) / \text{min} \times 8199$$

V_t: Reaction total volume 1.02 ml, V_s: sample volume 0.02 ml

ε: NADH molar absorptivity 6.22,

1000: transfer U/ml to U/L, L_p: Optical path length (cm)

REFERENCE RANGE

Sex	Conc. (U/L)
Male	24-195 U/L
Female	24-170 U/L

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Hemolysis specimen shall be avoided.
- The test is developed to determine creatine kinase concentrations up to 800U/L. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.

MeDiPro CREATINE KINASE TEST (CK) - UV kinetic, NAC method

- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.

- Waste management please refers to the local legal requirements.

- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)

- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[U/L]	SD [U/L]	CV[%]
Sample1	124	1.11	0.89
Sample2	417	2.12	0.51
Sample3	452	4.37	0.97

- Linearity

This kit has a good linearity up to 1600U/L.

- Interference

Interference	Influence effect
Hemoglobin	Not suitable when hemolysis occur
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 1.4%

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C, Part 2. Clin. Chem. Lab. Med. 2002, 40: 631.
- New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clinica Chimica. Acta. 2003, 327: 69.

PARAMETER SETUP

Hitachi 7170/917 Applications

TEST	[CK]
ASSAY CODE	[RateA]:[24]-[34]
SAMPLE VOLUME	[4]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH(nm)	[405][340]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[CK]
ASSAY CODE	[RateA]:[35]-[50]
SAMPLE VOLUME	[6]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH(nm)	[405][340]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0015M	MeDiPro CREATINE KINASE TEST	R1 6×20ml R2 3×10ml
BC-0015A	MeDiPro CREATINE KINASE TEST	R1 4×60ml R2 2×30ml
BC-0015B	MeDiPro CREATINE KINASE TEST	R1 4×100ml R2 2×50ml
BC-0015C	MeDiPro CREATINE KINASE TEST R1	R1 2×300ml
BC-0015D	MeDiPro CREATINE KINASE TEST R1	R1 2×500ml
BC-0015G	MeDiPro CREATINE KINASE TEST R2	R2 2×200ml

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan