台塑肌酐試劑（CREA）- Picric acid test

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**MeDiPro CREATININE TEST (CREA) - Picric acid test****INTENDED USE**

For the quantitative determination of creatinine in serum or urine.

CLINICAL SIGNIFICANCE

Creatinine is excreted as a waste product by the kidneys. Increased serum creatinine level usually indicates impairment of renal function. Creatinine appears in the glomerular filtrate and is not reabsorbed by the tubule. Hence, any condition that reduces the glomerular filtration rate will result in the increase of creatinine concentration in plasma. Since the excretion rate of creatinine is relatively constant and its production rate is not influenced by protein catabolism or other external factors. The concentration of creatinine in serum is, therefore, a good measure of renal glomerular function

PRINCIPLE

Creatinine + picric acid $\xrightarrow{\text{Alkali}}$ Creatinine-picric acid
(Reddish colored complex)

REAGENT

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R ₁ :	NaOH	2.67 mmol/L
R ₂ :	Picric acid	6.5 mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expire date. If stored at **2 – 8 °C.**, reagent should be protected from light and contamination should be avoided.

Do not freeze the reagent!

SPECIMEN COLLECTION AND PREPARATION

The choice of specimen is serum and the use of plasma is not recommended. Creatinine in serum is stable for 7 days at room temperature, 28 days in the refrigerator and at least 180 days when frozen. Hemolysis should be avoided.

For urine specimens, it should be diluted 1:100 with water. Avoid conditions which could allow bacterial growth. Add preservatives such as fluoride, thymol or boric acid during collection is recommended.

PROCEDURES

- Main wavelength : 505 nm
Sub. wavelength : 570nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains two reagents, ready to use.

	Control	Specimen	Blank
Specimen (ml)	---	0.1	---
Control (ml)	0.1	---	---
ddH ₂ O (ml)	---	---	0.1
R1 (ml)	0.8	0.8	0.8
Mix, 37°C incubate 5min			
R2 (ml)	0.2	0.2	0.2

Mix, incubate at 37°C for 30sec, and read the initial absorbance A₁ against reagent blank, then read end absorbance A₂ in every 30sec. A = A₅₀₅ - A₅₇₀.

CALCULATION

With standard or calibrator

$$\text{Creatinine (mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

REFERENCE RANGE

Age	Conc. (mg/dL)
Serum(male)	0.6-1.1
Serum(female)	0.5-0.9
Collected urine	0.5-1.5g/24hr

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition.
- The test is developed to determine creatinine concentrations up to 10mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2

**MeDiPro CREATININE TEST (CREA) - Picric acid test**

as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.


- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)
- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.
- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	1.70	0.00	0.00
Sample2	7.10	0.00	0.00
Sample3	7.08	0.04	0.60

- Linearity

This kit has a good linearity up to 20mg/dL.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 500mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 16mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 4mg/dL
Intrafat	No interference was observed by intrafat up to 2.0%

- Stability

Expire day	1 year
Open vial stability	14 day

REFERENCE

Narayanan S, Appelton H. Creatinine: a review. Clin Chem 1980; 26: 1119.

PARAMETER SETUPHitachi 7170/917 Applications

TEST	[CREA]
ASSAY CODE	[RateA]:[19]-[23]
SAMPLE VOLUME	[14]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH (nm)	[570][505]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[CREA]
ASSAY CODE	[RateA]: [30]-[37]
SAMPLE VOLUME	[18]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH(nm)	[570][505]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0017M	MeDiPro CREATININE TEST	R1 6x20ml R2 3x10ml
BC-0017A	MeDiPro CREATININE TEST	R1 4x60ml R2 2x30ml
BC-0017B	MeDiPro CREATININE TEST	R1 4x100ml R2 2x50ml
BC-0017C	MeDiPro CREATININE TEST R1	R1 2x300ml
BC-0017D	MeDiPro CREATININE TEST R1	R1 2x500ml
BC-0017G	MeDiPro CREATININE TEST R2	R2 2x200ml

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan