## 台塑葡萄糖試劑 ( GLU-HK ) - Enzymatic UV test without deproteinisation

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## 台塑葡萄糖试剂 ( GLU-HK ) - Enzymatic UV test without deproteinisation

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## MeDiPro GLUCOSE TEST-Hexokinase method (GLU-HK) - Enzymatic UV test without deproteinisation

### INTENDED USE



For the quantitative determination of glucose in serum or plasma.

### CLINICAL SIGNIFICANCE

The most common disease associated with abnormal carbohydrate metabolism is diabetes mellitus, with its accompanying high blood glucose levels. Other conditions which may also result in abnormal blood glucose levels include: disorders of the pituitary gland, hyperthyroidism, Cushing's disease, traumatic injury, convulsive disorders, mental stress and pheochromocytoma. Acute and chronic infection, eclampsia, hypertension and severe liver disease may also exhibit transitory elevation of blood glucose level. On the other hand, hyperinsulinism from either exogenous insulin overdose or from lesions of the pancreas can result in low level of blood glucose.

### PRINCIPLE

The reagent used here is based on the hexokinase (HK) – glucose-6-phosphate dehydrogenase (G6PDH) U.V. end point method. The reactions are as follows:

The increase in NADH concentration is directly proportional to the glucose concentration.

### REAGENT

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R <sub>1</sub> :	Buffer pH6.5	
	ATP	2 mmol/L
	NAD <sup>+</sup>	2 mmol/L
R <sub>2</sub> :	G6PDH	1500 U/L
	Hexokinase	1500 U/L

### STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at 2 – 8 °C., contamination should be avoided.

**Do not freeze the reagent!**

### SPECIMEN COLLECTION AND PREPARATION

Both serum and plasma samples can be used. For serum samples, collect whole blood and allow it to clot in clean test tube at room temperature. Separate and then transfer the serum without delay to a clean test tube. Do the test as soon as possible or store at 2~8°C to avoid degradation. For plasma specimens, collect whole blood into a tube containing a suitable anticoagulant, (EDTA, heparin, etc.), separate and transfer the plasma into a clear test tube. To prevent degradation from glycolysis, fluoride (up to 10 mg/dl) may be added with no effect on the test results.

### PROCEDURES

- Main wavelength : 340 nm  
Sub. wavelength : 400nm  
Reaction Temperature : 37°C  
Optical path length : 1.0 cm
- This kit contains two reagents, ready to use.

	Blank	Control	Sample
Sample(ml)	---	---	0.01
Control (ml)	---	0.01	---
ddH <sub>2</sub> O (ml)	0.01	---	---
R <sub>1</sub> (ml)	0.8	0.8	0.8
Mix, 37°C incubate 5min			
R <sub>2</sub> (ml)	0.2	0.2	0.2

Mix, incubate at 37°C for 1 min, and read the initial absorbance A<sub>1</sub> against reagent blank, then read end absorbance A<sub>2</sub> in every 1 min. A = A<sub>340</sub> - A<sub>405</sub>.

### CALCULATION

With standard or calibrator

$$\text{Glucose(mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

### REFERENCE RANGE

70-105 mg/dL (3.9-5.8 mmol/L)

### WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition.
- The test is developed to determine glucose concentrations up to 400mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.

## MeDiPro GLUCOSE TEST-Hexokinase method (GLU-HK) - Enzymatic UV test without deproteinisation


- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)
- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.
- Validity please see the reagent box label shown.

### REAGENT CHARACTERS

#### 1. Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	89	0.94	1.06
Sample2	281	1.64	0.58
Sample3	277	1.75	0.63

#### 2. Linearity

This kit has a good linearity up to 600mg/dL.

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 500mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 2.0%

#### 4. Stability

Expire day	1 year
Open vial stability	30 day

### REFERENCE

- Henry, J.B., "Clinical Diagnosis and Management by Laboratory Method." W.B. Saunders and Company Philadelphia, PA, p. 153 (1979).
- Barthelmai, W., and Czek, R., Klin. Wochenscht., 40:585 (1962).
- Tietz, N.W., Fundamentals of Clinical Chemistry, 2 nd. Ed., W.B. Saunders Co., Philadelphia, PA243 (1976).

### PARAMETER SETUP

#### Hitachi 7170 / 917 Applications

TEST	[GLU-HK]
ASSAY CODE	[2 POINT]: [16]-[34]
SAMPLE VOLUME	[2]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH (nm)	[405][340]
CALIB. METHOD	[Linear]

#### Hitachi 7150 / 717 Applications

TEST	[GLU-HK]
ASSAY CODE	[2 POINT]: [24]-[50]
SAMPLE VOLUME	[3]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH (nm)	[405][340]
CALIB. METHOD	[Linear]

### ORDERING INFORMATION

Cat. No.	Product	Package
BC-0032M	MeDiPro GLUCOSE TEST - Hexokinase Method	R1 6x20ml R2 3x10ml
BC-0032A	MeDiPro GLUCOSE TEST - Hexokinase Method	R1 4x60ml R2 2x30ml
BC-0032B	MeDiPro GLUCOSE TEST - Hexokinase Method	R1 4x100ml R2 2x50ml
BC-0032C	MeDiPro GLUCOSE TEST - Hexokinase Method R1	R1 4x300ml
BC-0032D	MeDiPro GLUCOSE TEST - Hexokinase Method R1	R1 4x500ml
BC-0032G	MeDiPro GLUCOSE TEST - Hexokinase Method R2	R2 4x200ml

### FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

### FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan