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
MeDiPro GLUCOSE TEST-Oxidase method(GLU-OD) - Enzymatic colorimetric test (GOD-PAP)


INTENDED USE

For the quantitative determination of glucose in serum or plasma.

CLINICAL SIGNIFICANCE

The most common disease associated with abnormal carbohydrate metabolism is diabetes mellitus, with its accompanying high blood glucose levels. Other conditions which may also result in abnormal blood glucose levels include: disorders of the pituitary gland, hyperthyroidism, Cushing's disease, traumatic injury, convulsive disorders, mental stress and pheochromocytoma. Acute and chronic infection, eclampsia, hypertension and severe liver disease may also exhibit transitory elevation of blood glucose level. On the other hand, hyperinsulinism from either exogenous insulin overdose or from lesions of the pancreas can result in low level of blood glucose.

PRINCIPLE

REAGENT

- Package: please see the reagent box label shown.
- Components:

Component	Conc.
Buffer pH7.0	
4-Aminoantipyrine	1mmol/L
phenol	10mmol/L
glucose oxidase	36000U/L
peroxidase	12000U/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at 2 – 8 °C., reagent should be protected from light and contamination should be avoided. **Do not freeze the reagent!**

SPECIMEN COLLECTION AND PREPARATION

Both serum and plasma samples can be used. For serum samples, collect whole blood and allow it to clot in clean test tube at room temperature. Separate and then transfer the serum without delay to a clean test tube. Do the test as soon as possible or store at 2~8 °C to avoid degradation. For plasma specimens, collect whole blood into a tube containing a suitable anticoagulant, (EDTA, heparin, etc.), separate and transfer the plasma into a clear test tube. To prevent degradation from glycolysis, fluoride (up to 10 mg/dl) may be added with no effect on the test results.

PROCEDURES

- Main wavelength : 505 nm
Sub. wavelength : 660nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains single reagent, ready to use.

	Blank	Control	Specimen
Specimen(ml)	---	---	0.01
Control(ml)	---	0.01	---
ddH ₂ O(ml)	0.01	---	---
R1(ml)	1.0	1.0	1.0

Mix, incubate at 37°C for 5 min, and read the absorbance against reagent blank. $A = A_{505} - A_{660}$

CALCULATION

With standard or calibrator

$$\text{Glucose (mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

REFERENCE RANGE

70-105 mg/dL (3.9-5.8 mmol/L)

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. **Bio-Rad Lyphochek control is recommended to use as serum control.**
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition.
- The test is developed to determine glucose concentrations up to 400mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

MeDiPro GLUCOSE TEST-Oxidase method(GLU-OD) - Enzymatic colorimetric test (GOD-PAP)

potentially hazardous components. (MSDS could be obtained from local dealer.)

- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	93	0.83	0.90
Sample2	284	2.54	0.89
Sample3	280	2.22	0.79

- Linearity

This kit has a good linearity up to 600mg/dL.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 500mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 20mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 28mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 20mg/dL
Intrafat	No interference was observed by intrafat up to 0.2%

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

- Bishop ML, Duben-Engelkirk JL, Fody EP. Clinical Chemistry: principles, procedures, correlations. Fourth edition, 2000.
- Tietz, N.W., Fundamentals of Clinical Chemistry, 2nd. Ed., W.B. Saunders Co., Philadelphia, PA243 (1976).

PARAMETER SETUP

Hitachi 7170 / 917 Applications

TEST	[GLU-OD]
ASSAY CODE	[1 POINT]: [34]-[0]
SAMPLE VOLUME	[2]
R1 VOLUME	[200]
R2 VOLUME	[0]
WAVELENGTH (nm)	[660][505]
CALIB. METHOD	[Linear]

Hitachi 7150 / 717 Applications

TEST	[GLU-OD]
ASSAY CODE	[1 POINT]: [50]-[0]
SAMPLE VOLUME	[3]
R1 VOLUME	[300]
R2 VOLUME	[0]
WAVELENGTH (nm)	[660][505]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0019M	MeDiPro GLUCOSE TEST - Oxidase method	R1 6x20ml
BC-0019A	MeDiPro GLUCOSE TEST - Oxidase method	R1 4x60ml
BC-0019B	MeDiPro GLUCOSE TEST - Oxidase method	R1 4x100ml
BC-0019C	MeDiPro GLUCOSE TEST - Oxidase method	R1 4x300ml
BC-0019D	MeDiPro GLUCOSE TEST - Oxidase method	R1 4x500ml

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan