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
**MeDiPro HDL-CHOLESTEROL TEST - Direct method****INTENDED USE**

For the quantitative determination of HDL-cholesterol in serum.

**CLINICAL SIGNIFICANCE**

HDL-cholesterol(HDL-C) is a very strong negative correlation between the risk of coronary artery disease. Low levels predict high risk and vice versa. Results from Framingham study showed higher incidence rates of coronary artery disease when HDL-cholesterol was below 45 mg/dl. However, HDL-C levels do not necessarily correlate with total cholesterol and LDL-cholesterol. A person may have normal total cholesterol content yet at high risk due to too low levels of HDL-cholesterol.

**PRINCIPLE**

CHOD: cholesterol oxidase  
CHOE: cholesterol esterase  
POD: peroxidase

**REAGENT**

3. Package: please see the reagent box label shown.  
4. Components:

	Component	Conc.
R <sub>1</sub> :	MOPS buffer	100mmol/L
	HDAOS	2mmol/L
	Magnesium sulfate	10mmol/L
R <sub>2</sub> :	PIPES	100mmol/L
	4-AA	0.5mmol/L
	CHOE	≥0.5KU/L
	CHOD	≥1KU/L
	POD	≥10KU/L

**STORE TEMPERATURE**

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C.**, reagent should be protected from light and contamination should be avoided.

**Do not freeze the reagent!**

**SPECIMEN COLLECTION AND PREPARATION**

Blood specimen should be drawn after the patient has fasted for at least 12 hours. Fresh serum and plasma are the choices.

**PROCEDURES**

1. Main wavelength : 600 nm  
Sub. wavelength : 700nm  
Reaction Temperature : 37°C  
Optical path length : 1.0 cm

2. This kit contains two reagents, ready to use.

	Specimen	Control	Blank
Specimen (ml)	0.004	---	---
Control (ml)	---	0.004	---
ddH <sub>2</sub> O (ml)	---	---	0.004
R <sub>1</sub> (ml)	0.3	0.3	0.3
Mix, 37 incubate 5min			
R <sub>2</sub> (ml)	0.1	0.1	0.1

Mix, incubate at 37°C for 5min, and read the absorbance against reagent blank at 600nm and 700nm.  $A = A_{600} - A_{700}$ .

**CALCULATION**

With standard or calibrator

$$\text{HDL-C (mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

**REFERENCE RANGE**

32~76 mg/dL

**WARNINGS AND PRECAUTIONS**

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Bilirubin might be interfere the result.
- The test is developed to determine HDL-C concentrations up to 150mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal

**MeDiPro HDL-CHOLESTEROL TEST - Direct method****requirements.**


9. Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)
10. According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.
11. Validity please see the reagent box label shown.

**REAGENT CHARACTERS**

1. Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	63.8	0.41	0.64
Sample2	36.4	0.27	0.73
Sample3	35.1	0.24	0.69

2. Linearity

This kit has a good linearity up to 250mg/dL.

3. Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 500mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 100mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 16mg/dL
Intrafat	No interference was observed by intrafat up to 3%

4. Stability

Expire day	1 year
Open vial stability	30 day

**REFERENCE**

1. Lopes-Virella, M.F., Stone, P., Ellis, S., Colwell, J.A., Cholesterol Determination in High-Density Lipoproteins Separated by Three Different Methods, Clinical Chem. 23/5, 882-884, 1977.

2. Demacker, P.N.M., Hifmans, A.G.M., Vos-Janssen, H.E., van't Laar, A., Jansen, A.P., A Study of the Use of Polyethylene Glycol in Estimating Cholesterol in High-Density Lipoprotein, Clinical Chem. 26/13,1775-1779, 1980.

**PARAMETER SETUP****Hitachi 7170/917 Applications**

TEST	[HDL]
ASSAY CODE	[2 Point]: [16]-[34]
SAMPLE VOLUME	[2]
R1 VOLUME	[180]
R2 VOLUME	[60]
WAVELENGTH(nm)	[700][600]
CALIB. METHOD	[Linear]

**Hitachi 7150/717 Applications**

TEST	[HDL]
ASSAY CODE	[2 Point]: [24]-[50]
SAMPLE VOLUME	[3]
R1 VOLUME	[240]
R2 VOLUME	[80]
WAVELENGTH(nm)	[700][ 600]
CALIB. METHOD	[Linear]

**ORDERING INFORMATION**

Cat. No.	Product	Package
BC-0020M	MeDiPro HDL-CHOLESTEROL TEST	R1 6x15ml R2 3x10ml
BC-0020A	MeDiPro HDL-CHOLESTEROL TEST	R1 4x45ml R2 2x30ml
BC-0020B	MeDiPro HDL-CHOLESTEROL TEST	R1 3x90ml R2 3x30ml
BC-0020C	MeDiPro HDL-CHOLESTEROL TEST R1	R1 2x300ml
BC-0020D	MeDiPro HDL-CHOLESTEROL TEST R1	R1 2x500ml
BC-0020G	MeDiPro HDL-CHOLESTEROL TEST R2	R2 2x100ml

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan