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**MeDiPro LACTATE DEHYDROGENASE TEST (LDH) - LDH-P method****INTENDED USE**

For the quantitative determination of lactate dehydrogenase activity in serum.

**CLINICAL SIGNIFICANCE**

Lactate dehydrogenase (LDH) is widely distributed in mammalian tissues, being rich in myocardium, kidney, liver and muscle. Determination of serum LDH activity is one of the most frequently performed assays in the diagnosis of myocardial and pulmonary infarction. Other conditions, such as metaloblastic anemia, extensive carcinomatosis, severe shock and hypoxia, granulocytic or acute anemia, hemolytic anemia, infectious mononucleosis, progressive muscular dystrophy, hepatitis, cirrhosis, obstructive jaundice, and in delirium tremens all will cause the increased activity of LDH.

**PRINCIPLE**

The rate of NADH formation is proportional to the activity of LDH

**REAGENT**

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R <sub>1</sub> :	Buffer	(pH8.95)±0.3
	Pyruvate	55 mmol/L
R <sub>2</sub> :	NADH	7.5 mmol/L

**STORE TEMPERATURE**

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C.**, reagent should be protected from light and contamination should be avoided.  
**Do not freeze the reagent!**

**SPECIMEN COLLECTION AND PREPARATION**

Serum is the choice for the assay but heparin and EDTA which will inhibit the LDH activity should not be used. Hemolysis should be avoided. Serum LDH is reported to be stable for 1 week at room temperature, and for 3 weeks at 2~8 °C. Specimen frozen should be avoided since the activities of LDH-4 and LDH-5 decrease significantly after frozen.

**PROCEDURES**

- Main wavelength : 340 nm  
Sub. wavelength : 405nm  
Reaction Temperature : 37°C  
Optical path length : 1.0 cm

- This kit contains two reagents, ready to use.

Volume (ml)	
Sample	0.02
R <sub>1</sub>	0.8
Mix, 37°C incubate 5min	
R <sub>2</sub>	0.2

Mix, incubate at 37°C for 1 min, and read the initial absorbance A<sub>1</sub> against reagent blank, then read end absorbance A<sub>2</sub> in every 1 min. A=A<sub>340</sub>-A<sub>405</sub>.

**CALCULATION**

$$\text{LDH (U/L)} = \frac{(A_2 - A_1) / \text{min} \times V_t \times 1000}{L_p \times \epsilon \times V_s}$$
$$= (A_2 - A_1) / \text{min} \times 8199$$

V<sub>t</sub>: Reaction total volume 1.02 ml, V<sub>s</sub>: sample volume 0.02 ml  
ε: NADH molar absorptivity 6.22,  
1000: transfer U/ml to U/L, L<sub>p</sub>: Optical path length (cm)

**REFERENCE RANGE**

200 ~ 480 U/L

**WARNINGS AND PRECAUTIONS**

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Cause lots of LDH in red cell, so hemolysis should be avoided.
- The test is developed to determine LDH concentrations up to 1000U/L. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- Do not use if reagents become turbid or initial blank OD less than 0.6.
- Fresh serum must be separated in 2 hours, Hg, sodium borate, and most anticoagulant, except heparin, might inhibit LDH activity.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.

**MeDiPro LACTATE DEHYDROGENASE TEST (LDH) - LDH-P method**

- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.

- Waste management please refers to the local legal requirements.

- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)

- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


- Validity please see the reagent box label shown.

**REAGENT CHARACTERS**

- Precision (Within run)

N=15	Mean[U/L]	SD [U/L]	CV[%]
Sample1	290	2.53	0.87
Sample2	845	2.06	0.24
Sample3	9656	6.04	0.63

- Linearity

This kit has a good linearity up to 2000U/L.

- Interference

Interference	Influence effect
Hemoglobin	Not suitable when hemolysis occur
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 1.6%

- Stability

Expire day	1 year
Open vial stability	30 day

**REFERENCE**

- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C, Part 3. Clin. Chem. Lab. Med. 2002, 40: 631.
- New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clinica Chimica. Acta. 2003, 327: 69.

**PARAMETER SETUP**

Hitachi 7170 / 917 Applications

TEST	[LDH]
ASSAY CODE	[Rate A]: [19]-[30]
SAMPLE VOLUME	[4]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH (nm)	[405][340]
CALIB. METHOD	[Linear]

Hitachi 7150 / 717 Applications

TEST	[LDH]
ASSAY CODE	[Rate A]: [30]-[45]
SAMPLE VOLUME	[6]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH (nm)	[405][340]
CALIB. METHOD	[Linear]

**ORDERING INFORMATION**

Cat. No.	Product	Package
BC-0022M	MeDiPro LACTATE	R1 6×20ml
	DEHYDROGENASE TEST	R2 3×10ml
BC-0022A	MeDiPro LACTATE	R1 4×60ml
	DEHYDROGENASE TEST	R2 2×30ml
BC-0022B	MeDiPro LACTATE	R1 4×100ml
	DEHYDROGENASE TEST	R2 2×50ml
BC-0022C	MeDiPro LACTATE	R1 2×300ml
	DEHYDROGENASE TEST R1	
BC-0022D	MeDiPro LACTATE	R1 2×500ml
	DEHYDROGENASE TEST R1	
BC-0022G	MeDiPro LACTATE	R2 2×200ml
	DEHYDROGENASE TEST R2	

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan