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

**MeDiPro LDL-CHOLESTEROL TEST - Direct method****INTENDED USE**

For the quantitative determination of LDL-cholesterol in serum.

**CLINICAL SIGNIFICANCE**

Of all the forms of cholesterol in the blood, the LDL cholesterol is considered the most important form in determining risk of heart disease. Elevated levels of LDL indicate risk for coronary artery disease and arteriosclerosis.

**PRINCIPLE**

CHOD: cholesterol oxidase  
CHOE: cholesterol esterase  
POD: peroxidase

**REAGENT**

3. Package: please see the reagent box label shown.

4. Components:

Component	Conc.
R <sub>1</sub> : HDAOS	3mmol/L
R <sub>2</sub> : CHOD	≥ 0.5KU/L
CHOE	≥ 0.5KU/L
POD	≥ 10KU/L
4-AA	0.5mmol/L

**STORE TEMPERATURE**

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C**., reagent should be protected from light and contamination should be avoided.  
**Do not freeze the reagent!**

**SPECIMEN COLLECTION AND PREPARATION**

Blood specimen should be drawn after the patient has fasted for at least 12 hours. Fresh serum and plasma are the choices.

**PROCEDURES**

1. Main wavelength : 600 nm  
Sub. wavelength : 700nm  
Reaction Temperature : 37°C  
Optical path length : 1.0 cm

2. This kit contains two reagents, ready to use.

	Specimen	Control	Blank
Specimen (ml)	0.0025	---	---
Control (ml)	---	0.0025	---
ddH <sub>2</sub> O (ml)	---	---	0.0025
R <sub>1</sub> (ml)	0.24	0.24	0.24
Mix, 37 incubate 5min			
R <sub>2</sub> (ml)	0.08	0.08	0.08

Mix, incubate at 37°C for 5min, and read the absorbance against reagent blank at 600nm and 700nm.  $A = A_{600} - A_{700}$ .

**CALCULATION**

With standard or calibrator

$$\text{LDL-C (mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

**REFERENCE RANGE**

65~160mg/dL

**WARNINGS AND PRECAUTIONS**

1. This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
2. Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
3. This kit is for professionals and *in vitro* diagnostic use only.
4. To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Bilirubin might be interfere the result.
5. The test is developed to determine LDL-C concentrations up to 450mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
6. The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
7. Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.

**MeDiPro LDL-CHOLESTEROL TEST - Direct method**

8. Waste management please refers to the local legal requirements.

9. Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)

10. According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


11. Validity please see the reagent box label shown.

**REAGENT CHARACTERS**

1. Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	188.0	1.04	0.56
Sample2	90.9	1.09	1.19
Sample3	89.7	1.34	1.49

2. Linearity

This kit has a good linearity up to 500mg/dL.

3. Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 500mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 100mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 3%

4. Stability

Expire day	1 year
Open vial stability	30 day

**REFERENCE**

Bachorik PS. Measurement of low-density lipoprotein cholesterol. In: Rifai N, Warnick GR, eds. Handbook of lipoprotein testing. Washington, DC: AACC press, 1997; 145.

**PARAMETER SETUP**

Hitachi 7170/917 Applications

TEST	[LDL]
ASSAY CODE	[2 POINT]: [16]-[34]
SAMPLE VOLUME	[2]
R1 VOLUME	[180]
R2 VOLUME	[60]
WAVELENGTH(nm)	[700][600]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[LDL]
ASSAY CODE	[2 POINT]: [24]-[50]
SAMPLE VOLUME	[3]
R1 VOLUME	[240]
R2 VOLUME	[80]
WAVELENGTH(nm)	[700][600]
CALIB. METHOD	[Linear]

**ORDERING INFORMATION**

Cat. No.	Product	Package
BC-0023M	MeDiPro LDL-CHOLESTEROL TEST	R1 6x15ml R2 3x10ml
BC-0023A	MeDiPro LDL-CHOLESTEROL TEST	R1 4x45ml R2 2x30ml
BC-0023B	MeDiPro LDL-CHOLESTEROL TEST	R1 3x90ml R2 3x30ml
BC-0023C	MeDiPro LDL-CHOLESTEROL TEST R1	R1 2x300ml
BC-0023D	MeDiPro LDL-CHOLESTEROL TEST R1	R1 2x500ml
BC-0023G	MeDiPro LDL-CHOLESTEROL TEST R2	R2 2x100ml

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan