台塑鎂試劑(MG) - Calmagite test

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MeDiPro MAGNESIUM TEST (MG) - Calmagite test

INTENDED USE

For the quantitative determination of magnesium in serum.

CLINICAL SIGNIFICANCE

Magnesium is one of the most abundant cations in the body and is involved in many biochemical reactions. Many enzymes such as alkaline phosphatase, ALP require magnesium as an activator. Magnesium is also necessary for the stability of conformational structure of many macromolecules such as DNA, RNA, etc...

Although little is known about the regulation of magnesium levels in blood, it has been reported that para-thyroid gland is involved. Increased level of magnesium has been shown in Addison's disease, diabetic acidosis, renal failure and vitamin D intoxication, and decreased level of magnesium are observed in diabetes, diuretics, hyperthyroidism, hyperalimentation, alcoholism, myocardial infarction, congestive heart failure and liver cirrhosis.

PRINCIPLE

Magnesium forms a colored complex with calmagite under alkaline conditions with an absorbance at 505 nm. The interference of calcium is eliminated by the addition of EGTA. The reaction equation is shown as following:

REAGENT

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R ₁ :	EGTA	0.075mmol/L
	KCl	134 mmol/L
	KOH	176 mmol/L
R ₂ :	calmagite	0.064mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C.**, reagent should be protected from light and contamination should be avoided.
Do not freeze the reagent!

SPECIMEN COLLECTION AND PREPARATION

Freshly drawn serum is the specimen of choice, plasma derived from specimen collection tubes containing EDTA, citrate, or oxalate must not be used in this test. Serum should be removed from the clot without delay to avoid loss of serum magnesium due to increased erythrocyte permeability to magnesium.

PROCEDURES

- Main wavelength : 505 nm
Sub. wavelength : 700nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains two reagents, ready to use.

Volume (ml)	
Sample	0.01
R ₁	0.8
Mix, 37°C incubate 3min	
R ₂	0.2

Mix, incubate at 37°C for 1 min, and read the initial absorbance A₁ against reagent blank, then read end absorbance A₂ in every 1 min. A=A₅₀₅-A₇₀₀

CALCULATION

With standard or calibrator

$$\text{MG(mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std.}/\text{cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

REFERENCE RANGE

1.6~2.6 mg/dL (0.65~1.05 mmol/L)

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Lipemia specimen should be avoided.
- The test is developed to determine magnesium concentrations up to 6.0mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- Do not use if reagents become turbid or initial blank OD over 1.0.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.

MeDiPro MAGNESIUM TEST (MG) - Calmagite test

8. Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.

9. Waste management please refers to the local legal requirements.

10. Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)


11. According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.

12. Validity please see the reagent box label shown.

REAGENT CHARACTERS

1. Precision (Within run)			
N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	2.61	0.10	3.64
Sample2	4.63	0.04	0.94
Sample3	4.47	0.02	0.46

2. Linearity

This kit has a good linearity up to 8.0mg/dL.

3. Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 150mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	Not suitable when lipemia occur

4. Stability

Expire day	1 year
Open vial stability	14 day

REFERENCE

- Liedtke RJ., Kroon G. Automated calmagite compleximetric measurement of magnesium in serum, with sequential addition of EDTA to eliminate endogenous interference. Clin. Chem. 1984. 30:1801.
- Polancic JE. Magnesium: metabolism, clinical importance, and analysis. Clin Lab Sci 1991; 4(2).

PARAMETER SETUP

Hitachi 7170/917 Applications

TEST	[MG]
ASSAY CODE	[2 point end]:[16]-[34]
SAMPLE VOLUME	[2]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH(nm)	[700][505]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[MG]
ASSAY CODE	[2 point end]:[24]-[50]
SAMPLE VOLUME	[3]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH(nm)	[700][505]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0024M	MeDiPro MAGNESIUM TEST	R1 6×20ml R2 3×10ml
BC-0024A	MeDiPro MAGNESIUM TEST	R1 4×60ml R2 2×30ml
BC-0024B	MeDiPro MAGNESIUM TEST	R1 4×100ml R2 2×50ml
BC-0024C	MeDiPro MAGNESIUM TEST R1	R1 2×300ml
BC-0024D	MeDiPro MAGNESIUM TEST R1	R1 2×500ml
BC-0024G	MeDiPro MAGNESIUM TEST R2	R2 2×200ml

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan