## 台塑磷試劑 (PHOS) - Phosphomolybdate method

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## MeDiPro PHOSPHORUS TEST (PHOS) - Phosphomolybdate method


### INTENDED USE

For the quantitative determination of inorganic phosphorus in serum or plasma.

### CLINICAL SIGNIFICANCE

Phosphorus is an abundant element in the body and plays an important role in the intermediary metabolism, skeletal formation, dentition and acid-base balance.

### PRINCIPLE

### REAGENT

- Package: please see the reagent box label shown.
- Components:

Component	Conc.
Sulfuric acid	10 mmol/L
Ammonium molybdate	3.4 mmol/L

### STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at 2 – 8 °C., contamination should be avoided.

**Do not freeze the reagent!**

### SPECIMEN COLLECTION AND PREPARATION

Fasting blood should be used since glucose ingestion will lower the inorganic phosphorus level. Fresh serum and heparin-treated plasma are the choices. Fresh unhemolyzed serum should be separated as soon as possible from the red blood cells. Hemolyzed serum should be avoided since it contains organic phosphorus compounds which may be hydrolyzed and interferes with inorganic phosphorus. If the serum is properly separated, the phosphorus is stable for 1 week when stored refrigerated.

### PROCEDURES

- Main wavelength : 340 nm  
Sub. wavelength : 405nm  
Reaction Temperature : 37°C  
Optical path length : 1.0 cm
- This kit contains single reagent, ready to use.

	Blank	Control	Specimen
ddH <sub>2</sub> O	0.01	---	---
Control	---	0.01	---
Specimen	---	---	0.01
R1	1.0	1.0	1.0

Mix, incubate at 37°C for 5 min, and read the absorbance against reagent blank.  $A = A_{340} - A_{405}$

### CALCULATION

With standard or calibrator

$$\text{Phosphorus (mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

### REFERENCE RANGE

Age	Conc. (mg/dL)
Child	4.0~5.5 mg/dL
Adult	2.5~4.5 mg/dL

### WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Hemolysis, ascorbic acid, and bilirubin might interfere with the result.
- Each instrument, material and distilled water should avoid calcium contamination.
- The test is developed to determine phosphorus concentrations up to 13mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)

## MeDiPro PHOSPHORUS TEST (PHOS) - Phosphomolybdate method


- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.
- Validity please see the reagent box label shown.

### REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	4.22	0.09	2.18
Sample2	8.14	0.08	0.96
Sample3	8.88	0.09	0.97

- Linearity

This kit has a good linearity up to 15mg/dL.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 100mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 30mg/dL
Bilirubin (free form)	Not suitable when jaundice occur
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	Not suitable when lipemia occur

- Stability

Expire day	1 year
Open vial stability	14 day

### REFERENCE

- Concustell E, Cortes M, Ferragut A, Gener J. Inorganic phosphorus measurement: an improved method. Clin Chim Acta. 1977 Dec 15;81(3):267-72.
- Bishop M.L., Duben-Engelkirk J.L., Fody E.P. Clinical Chemistry Fourth edition.

### PARAMETER SETUP

#### Hitachi 7170/917 Applications

TEST	[PHOS]
ASSAY CODE	[1 POINT]:[34]-[0]
SAMPLE VOLUME	[2]
R1 VOLUME	[200]
R2 VOLUME	[0]
WAVELENGTH(nm)	[405][340]
CALIB. METHOD	[Linear]

#### Hitachi 7150/717 Applications

TEST	[PHOS]
ASSAY CODE	[1 POINT]:[50]-[0]
SAMPLE VOLUME	[3]
R1 VOLUME	[300]
R2 VOLUME	[0]
WAVELENGTH(nm)	[405][340]
CALIB. METHOD	[Linear]

### ORDERING INFORMATION

Cat. No.	Product	Package
BC-0025M	MeDiPro PHOSPHORUS TEST	R1 6x20ml
BC-0025A	MeDiPro PHOSPHORUS TEST	R1 4x60ml
BC-0025B	MeDiPro PHOSPHORUS TEST	R1 4x100ml
BC-0025C	MeDiPro PHOSPHORUS TEST	R1 3x250ml
BC-0025D	MeDiPro PHOSPHORUS TEST	R1 2x500ml

### FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

### FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan