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## 台塑總膽紅素試劑（TBIL）- Diazo-DPD method

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
For the quantitative determination of total bilirubin in serum or plasma.

### CLINICAL SIGNIFICANCE

Bilirubin, is derived primarily in bile, from the breakdown of hemoglobin when senescent red blood cells are phagocytized. Normally, about 6 to 6.5 g of hemoglobin in aged red blood cells are broken down daily in an adult to form about 220 mg of bilirubin; another 50 to 60 mg of bilirubin originates from other sources. The detail mechanism of bilirubin formation is not well understood. But it is known that the heme group of hemoglobin is converted to bilirubin in reticuloendothelial system, which binds to albumin in plasma and is esterified to bilirubin diglucuronide (BDG), which is secreted from liver as a waste product. It presents in serum in the free and conjugated forms. An increase in the formation or retention of bilirubin in the body results in the increased levels of serum bilirubin and jaundice. This hyperbilirubinemia is classified as prehepatic, hepatic or posthepatic depending on the major cause of condition. Therefore, the determination of total bilirubin and its conjugated (direct) bilirubin is important for the diagnosis of hyperbilirubinemia.

### PRINCIPLE

### REAGENT

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R <sub>1</sub>	Sulfamic acid	70 mmol/L
	HCl	70 mmol/L
	DMSO	5 mol/L
R <sub>2</sub>	3,5-DCPD	5mmol/L

### STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C**., reagent should be protected from light and contamination is avoided.

**Do not freeze the reagent!**

### SPECIMEN COLLECTION AND PREPARATION

Fresh serum and plasma are the choices. The sample be collected without hemolysis, since hemoglobin inhibits the diazo reaction.

Avoid direct light exposure to the specimen since Bilirubin values may decrease as much as 50 % in one hour. Serum specimen may be kept in dark in 2-8°C for up to one week, and in freezer for 3 months without appreciable change in the Bilirubin levels.

### PROCEDURES

- Main wavelength : 546 nm  
Sub. wavelength : 660nm  
Reaction Temperature : 37°C  
Optical path length : 1.0 cm

- This kit contains two reagents, ready to use.

	Blank	Control	Specimen
R <sub>1</sub> (ml)	0.8	0.8	0.8
ddH <sub>2</sub> O(ml)	0.02	---	---
Control (ml)	---	0.02	---
Specimen(ml)	---	---	0.02

mix, 37°C incubate 5min

R <sub>2</sub> (ml)	0.2	0.2	0.2
---------------------	-----	-----	-----

Mix, incubate at 37°C for 5 min, and read the absorbance against reagent blank.  $A = A_{546} - A_{660}$

### CALCULATION

With standard or calibrator

$$\text{Total bilirubin} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

### REFERENCE RANGE

Age	Conc. (mg/dL)
Newborn	< 12
Adult	< 1.2

### WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition.
- Do not use if reagent 2 turns yellow.
- The test is developed to determine bilirubin concentrations up to 18mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2

## MeDiPro TOTAL BILIRUBIN TEST ( TBIL ) - Diazo-DPD method

as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.

- Waste management please refers to the local legal requirements.

- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)

- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


- Validity please see the reagent box label shown.

### REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	0.74	0.01	1.78%
Sample2	4.89	0.10	1.97%
Sample3	4.77	0.08	1.67%

- Linearity

This kit has a good linearity up to 20mg/dL.

- Interference

Interference	Influence effect
Hemoglobin	Not suitable when hemolysis occur
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Intrafat	Not suitable when lipemia occur

- Stability

Expire day	1 year
Open vial stability	30 day

### REFERENCE

- Basll T.D., Billy W.P., Edward A.S., Jon V.S. Standardization in bilirubin assays: evaluation of selected methods and stability of bilirubin solutions. 1973, 19: 984-993.
- Billy W.P., Basll T.D., etc. A candidate reference method for determination of bilirubin in serum. Test for transferability. Clin. Chem. 1983, 29: 297-301.

### PARAMETER SETUP

#### Hitachi 7170/917 Applications

TEST	[TBIL]
ASSAY CODE	[2POINT]: [16]-[34]
SAMPLE VOLUME	[4]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH(nm)	[660][546]
CALIB. METHOD	[Linear]

#### Hitachi 7150/717 Applications

TEST	[TBIL]
ASSAY CODE	[2POINT]: [24]-[50]
SAMPLE VOLUME	[6]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH(nm)	[660][546]
CALIB. METHOD	[Linear]

### ORDERING INFORMATION

Cat. No.	Product	Package
BC-0010M	MeDiPro TOTAL BILIRUBIN TEST	R1 6x20ml R2 3x10ml
BC-0010A	MeDiPro TOTAL BILIRUBIN TEST	R1 4x60ml R2 2x30ml
BC-0010B	MeDiPro TOTAL BILIRUBIN TEST	R1 4x100ml R2 2x50ml
BC-0010C	MeDiPro TOTAL BILIRUBIN TEST	R1 2x300ml
BC-0010D	MeDiPro TOTAL BILIRUBIN TEST	R1 2x500ml
BC-0010G	MeDiPro TOTAL BILIRUBIN TEST	R1 R2 2x200ml

### FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

### FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan