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

**MeDiPro TRIGLYCERIDE TEST (TG) - Enzymatic method****INTENDED USE**

For the quantitative determination of triglycerides in serum or plasma.

CLINICAL SIGNIFICANCE

Triglycerides (TG) circulate in blood associated with proteins in the form of lipoproteins. Elevated level of triglycerides in plasma was reported as a possible indicator of atherosclerosis. Diseases such as nephrotic syndrome, glycogen storage disease, biliary cirrhosis and others resulting from derangements in lipid metabolism also show increased triglyceride levels in serum.

PRINCIPLE**REAGENT**

- Package: please see the reagent box label shown.
- Components:

Component	Conc.
Buffer pH6.5±0.25	30 mmol/L
ATP	3 mmol/L
4-chlorophenol	0.8 mmol/L
NaOH	20 mmol/L
Glycerol kinase	400 U/L
4-AA	0.6 mmol/L
GPOD	7000 U/L
Peroxidase	1500 U/L
Lipase	200000 U/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C**., reagent should be protected from light and contamination should be avoided.

Do not freeze the reagent!

SPECIMEN COLLECTION AND PREPARATION

Fresh serum and plasma are the choices. Blood specimen should be drawn after the patient has fasted for at least 12 hours. All tubes and collection materials should be free of glycerol. Prolonged storage of samples at room temperature should be avoided in order to minimize possible degradation of triglyceride by endogenous lipase or by lipase secreted from contaminated bacteria.

PROCEDURES

- Main wavelength : 505 nm
Sub. wavelength : 660nm
Reaction Temperature : 37°C

Optical path length : 1.0 cm			
2. This kit contains single reagent, ready to use.			
	Blank	Control	Specimen
Specimen (ml)	—	—	0.01
Control (ml)	—	0.01	—
ddH ₂ O (ml)	0.01	—	—
R1 (ml)	1.0	1.0	1.0

Mix, incubate at 37°C for 10 min, and read the absorbance against reagent blank. $A = A_{505} - A_{660}$

CALCULATION

With standard or calibrator

$$\text{TG(mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

REFERENCE RANGE

36-165 mg/dL (0.41-1.86 mmol/L)

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Bilirubin and lipemia might be interfere the result.
- The test is developed to determine triglycerides concentrations up to 1800U/L. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal requirements.

**MeDiPro TRIGLYCERIDE TEST (TG) - Enzymatic method**

- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)

- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	201	1.62	0.80
Sample2	87	1.19	1.36
Sample3	94	1.01	1.08

- Linearity

This kit has a good linearity up to 2000mg/dL.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 50mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 40mg/dL
Bilirubin (free form)	Not suitable when jaundice occur
Bilirubin (conjugate form)	Not suitable when jaundice occur

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

- Klotzsch SG, McNamara JR. Triglyceride measurements: a review of methods and interferences. Clin. Chem. 1990. 36: 1605.
- McGowan M, Artiss J, Strandbergh DR, et al. A peroxidase-coupled method for the colorimetric determination of serum triglyceride. Clin. Chem. 1983. 29: 538

PARAMETER SETUP**Hitachi 7170/917 Applications**

TEST	[TG]
ASSAY CODE	[1 POINT]: [34]-[0]
SAMPLE VOLUME	[2]
R1 VOLUME	[200]
R2 VOLUME	[0]
WAVELENGTH (nm)	[660][505]
CALIB. METHOD	[Linear]

Hitachi 7150 / 717 Applications

TEST	[TG]
ASSAY CODE	[1 POINT]: [50]-[0]
SAMPLE VOLUME	[3]
R1 VOLUME	[300]
R2 VOLUME	[0]
WAVELENGTH (nm)	[660][505]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0027M	MeDiPro TRIGLYCERIDE TEST	R1 6×20ml
BC-0027A	MeDiPro TRIGLYCERIDE TEST	R1 4×60ml
BC-0027B	MeDiPro TRIGLYCERIDE TEST	R1 4×100ml
BC-0027C	MeDiPro TRIGLYCERIDE TEST	R1 4×300ml
BC-0027D	MeDiPro TRIGLYCERIDE TEST	R1 4×500ml

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

**FORMOSA BIOMEDICAL TECHNOLOGY CORP.**

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan