台塑總蛋白試劑（TP）- Biuret method

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MeDiPro TOTAL PROTEIN TEST (TP) - Biuret method


INTENDED USE

For the quantitative determination of Total Protein in serum or plasma.

CLINICAL SIGNIFICANCE

Proteins form the major portion of the solutes dissolved in the plasma fluid. The determination of total protein is useful in detecting hyperproteinemia due to homoconcentration will occur in dehydration, paraproteinemia, or monoclonal disease (multiple myeloma, macroglobulinemia, cryoglobulinemia) and in some chronic polyclonal diseases, liver cirrhosis, sarcoidosis, systemic lupus erythematosus, and chronic infections. Conditions will result in serum protein decrease including overhydration, protein loss through kidneys, severe burns, or in failure of protein synthesis (starvation, protein malnutrition, liver cell damage). Serum protein determination is also useful when determining the calcium concentration because the nondiffusible calcium fraction is bound to protein and varies with the serum protein.

PRINCIPLE

REAGENT

- Package: please see the reagent box label shown.
- Components:

Component	Conc.
CuSO ₄	12 mmol/L
NaOH	0.6 mmol/L
Potassium sodium tartrate	32 mmol/L
KI	30 mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C**., reagent should be protected from light and contamination should be avoided.

Do not freeze the reagent!

SPECIMEN COLLECTION AND PREPARATION

Serum and plasma are the specimen of choices. Avoid from hemolysis. Lipemic sera may cause elevated results. Sera with BSP will cause falsely elevated results. Proteins in serum are stable for 1 week at room temperature, up to one month in the refrigerator and longer when stored in the freezer.

PROCEDURES

- Main wavelength : 546 nm
Sub. wavelength : 700nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains single reagent, ready to use.

	Blank	Control	Specimen
ddH ₂ O (ml)	0.02	---	---
Control (ml)	---	0.02	---
Specimen (ml)	---	---	0.02
R1 (ml)	1.0	1.0	1.0

Mix, incubate at 37°C for 10 min, and read the absorbance against reagent blank. $A = A_{546} - A_{700}$

CALCULATION

With standard or calibrator

$$\text{Protein (g/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (g/dL)}$$

REFERENCE RANGE

6.0-8.2 (g/dL)

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Hemolysis and lipemia might interfere the test result.
- The test is developed to determine total protein concentrations up to 10 g/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)
- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

MeDiPro TOTAL PROTEIN TEST (TP) - Biuret method

the result.


- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[g/dL]	SD [g/dL]	CV[%]
Sample1	4.11	0.03	0.77
Sample2	4.16	0.03	0.77
Sample3	6.54	0.05	0.73

- Linearity

This kit has a good linearity up to 12g/dL.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 50mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 12mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 0.2%

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

Doumas BT, Bayse D, Borner K, et al. A candidate reference method for determination of total protein in serum : I. development and validation. Clin. Chem. 1981. 27:1642.

PARAMETER SETUP

Hitachi 7170/917 Applications

TEST	[TP]
ASSAY CODE	[1 POINT]:[34]-[0]
SAMPLE VOLUME	[4]
R1 VOLUME	[200]
R2 VOLUME	[0]
WAVELENGTH(nm)	[700][546]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[TP]
ASSAY CODE	[1 POINT]:[50]-[0]
SAMPLE VOLUME	[6]
R1 VOLUME	[300]
R2 VOLUME	[0]
WAVELENGTH(nm)	[700][546]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0026M	MeDiPro TOTAL PROTEIN TEST	R1 6x20ml
BC-0026A	MeDiPro TOTAL PROTEIN TEST	R1 4x60ml
BC-0026B	MeDiPro TOTAL PROTEIN TEST	R1 4x100ml
BC-0026C	MeDiPro TOTAL PROTEIN TEST	R1 4x300ml
BC-0026D	MeDiPro TOTAL PROTEIN TEST	R1 4x500ml

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan