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台塑尿酸試劑 (UA) - Enzymatic colorimetric test (uricase-PAP)

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
MeDiPro URIC ACID TEST (UA) - Enzymatic colorimetric test (uricase-PAP)

INTENDED USE
For the quantitative determination of uric acid in serum or plasma.

CLINICAL SIGNIFICANCE
Uric acid is the end-product of urine metabolism in human. It circulates and is filtered from the blood by the glomeruli. The determination of serum uric acid for detecting hyperuricemia is helpful in the diagnosis of gout, increased metabolism of nucleoproteins, such as in leukemia and polycythemia. Serum uric acid levels also increase with the decrease of renal function.

REAGENT

1. Package: please see the reagent box label shown.
2. Components:

Component	Conc.
R ₁ DHBS	0.01 mmol/L
Phosphate buffer pH6.7±0.25	150 mmol/L
R ₂ Uricase	200 U/L
Peroxidase	1000 U/L
4-aminoantipyrine	0.3 mmol/L

STORE TEMPERATURE
The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C.**, reagent should be protected from light and contamination should be avoided. **Do not freeze the reagent!**

SPECIMEN COLLECTION AND PREPARATION
Freshly drawn un-hemolyzed serum, plasma are the choice of specimen. Serum uric acid is stable 3 days at room temperature, one week when refrigerated; 6 months when frozen. Avoid bacteria contamination.

- PROCEDURES**
1. Main wavelength : 546 nm
Sub. wavelength : 700nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
 2. This kit contains two reagents, ready to use.

	Control	Specimen	Control
Specimen(ml)	—	0.02	—
Control(ml)	0.02	—	—
ddH ₂ O(ml)	—	—	0.02
R1 (ml)	0.8	0.8	0.8
Mix, 37 incubate 5min			
R2 (ml)	0.2	0.2	0.2

Mix, incubate at 37°C for 5 min, and read the absorbance against reagent blank. A=A₅₄₆-A₇₀₀

CALCULATION
With standard or calibrator

Uric acid (mg/dL)= $\frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$

REFERENCE RANGE
1.5-7.0 mg/dL (90-420 μmol/L)

- WARNINGS AND PRECAUTIONS**
1. This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
 2. Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
 3. This kit is for professionals and *in vitro* diagnostic use only.
 4. To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Bilirubin might interfere with the test result.
 5. The test is developed to determine uric acid concentrations up to 30mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
 6. The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
 7. Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.

FORMOSA BIOMEDICAL TECHNOLOGY CORP.
F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan


MeDiPro URIC ACID TEST (UA) - Enzymatic colorimetric test (uricase-PAP)

8. Waste management please refers to the local legal requirements.
9. Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)
10. According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.
11. Validity please see the reagent box label shown.

REAGENT CHARACTERS

1. Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	5.21	0.06	1.09
Sample2	9.99	0.09	0.88
Sample3	9.90	0.09	0.95

3. Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 500mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 40mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 24mg/dL
Bilirubin (conjugate form)	Not suitable when jaundice occur
Intrafat	No interference was observed by intrafat up to 2.0%

4. Stability

Expire day	1 year
Open vial stability	30 day

- REFERENCE**
- 1 · Kabasakalian P, Kalliney S, Wescott A. Determination of uric acid in serum, with the use of uricase and a tribromophenol- aminoantipyrine chromogen. Clin. Chem. 1973. 19: 522.
 - 2 · Gochman N, Schmitz JM. Automated determination of uric acid with use of a uricase-peroxidase system. Clin. Chem. 1971. 17: 1154.

PARAMETER SETUP
Hitachi 7170/917 Applications

TEST	[UA]
ASSAY CODE	[2 POINT]: [16]-[34]
SAMPLE VOLUME	[4]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH(nm)	[700][546]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[UA]
ASSAY CODE	[2 POINT]: [24]-[50]
SAMPLE VOLUME	[6]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH(nm)	[700][546]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0028M	MeDiPro URIC ACID TEST	R1 6×20ml R2 3×10ml
BC-0028A	MeDiPro URIC ACID TEST	R1 4×60ml R2 2×30ml
BC-0028B	MeDiPro URIC ACID TEST	R1 4×100ml R2 2×50ml
BC-0028C	MeDiPro URIC ACID TEST R1	R1 4×300ml
BC-0028D	MeDiPro URIC ACID TEST R1	R1 4×500ml
BC-0028G	MeDiPro URIC ACID TEST R2	R2 4×200ml

FORMOSA BIOMEDICAL TECHNOLOGY CORP.
F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan