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MeDiPro ALANINE AMINOTRANSFERASE (ALT) - UV-IFCC method


INTENDED USE

For the quantitative determination of alanine aminotransferase activity in serum.

CLINICAL SIGNIFICANCE

Alanine aminotransferase (ALT) is widely distributed in animal tissue, with the highest concentration present in liver and kidney. Acute destruction of tissue results in the release of ALT into the blood. Serum ALT level up to 1000 times normal has been observed in cases of acute hepatitis. Less elevated values have been associated with infectious mononucleosis, extra-hepatic, obstructive jaundice and cirrhosis.

PRINCIPLE

REAGENT

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R ₁	Tris Buffer	100 mmol/L
	LDH	>1200 U/L
	L-alanine	500 mmol/L
R ₂	α-ketoglutaric acid	15 mmol/L
	Tris Buffer	100 mmol/L
	NADH	0.18 mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at 2 – 8 °C., reagent should be protected from light and contamination should be avoided. **Do not freeze the reagent!**

SPECIMEN COLLECTION AND PREPARATION

Fresh serum is suitable for ALT assay. Avoid hemolysis. Hemolysis could cause false elevation of the ALT activity. Serum should be removed from the clot or cells without delay. Serum ALT activity is not stable and should be kept at 4°C for short period or stored frozen at -20°C or lower for longer periods.

PROCEDURES

- Main wavelength : 340 nm
Sub. wavelength : 405nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains two reagents, ready to use.

	Volume (ml)
sample	0.05
R ₁	0.8
Mix, 37°C incubate 5min	
R ₂	0.2

Mix, incubate at 37°C for 1 min, and read the initial absorbance A₁ against reagent blank, then read end absorbance A₂ in every 1 min. A=A₃₄₀-A₄₀₅.

CALCULATION

$$\begin{aligned} \text{ALT (U/L)} &= \frac{(A_2 - A_1) / \text{min} \times V_t \times 1000}{L_p \times \epsilon \times V_s} \\ &= (A_2 - A_1) / \text{min} \times 3376 \end{aligned}$$

V_t: Reaction total volume 1.05 ml, V_s: sample volume 0.05 ml
ε: NADH molar absorptivity 6.22,
1000: transfer U/ml to U/L, L_p: Optical path length (cm)

REFERENCE RANGE

<40U/L

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition.
- Do not use when reagent blank OD less than 1.0.
- The test is developed to determine alanine aminotransferase concentrations up to 800U/L. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

MeDiPro ALANINE AMINOTRANSFERASE (ALT) - UV-IFCC method

and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)

- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.


- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[U/L]	SD [U/L]	CV[%]
Sample1	32.90	0.74	2.24
Sample2	97.40	0.70	0.72
Sample3	97.40	0.70	0.72

- Linearity

This kit has a good linearity up to 1000U/L.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 500mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 30mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 12mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 2.0%

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C, Part 4. Clin. Chem. Lab. Med. 2002, 40: 631.
- New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clinica Chimica. Acta. 2003, 327: 69.

PARAMETER SETUP

Hitachi 7170/917 Applications

TEST	[ALT]
ASSAY CODE	[Rate A]: [19]-[34]
SAMPLE VOLUME	[10]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH(nm)	[405][340]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[ALT]
ASSAY CODE	[Rate A]: [30]-[50]
SAMPLE VOLUME	[15]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH(nm)	[405][340]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0007M	MeDiPro ALANINE AMINOTRANSFERASE TEST	R1 6×20ml, R2 3×10ml
BC-0007A	MeDiPro ALANINE AMINOTRANSFERASE TEST	R1 4×60ml, R2 2×30ml
BC-0007B	MeDiPro ALANINE AMINOTRANSFERASE TEST	R1 4×100ml, R2 2×50ml
BC-0007C	MeDiPro ALANINE AMINOTRANSFERASE TEST R1	R1 4×300ml
BC-0007D	MeDiPro ALANINE AMINOTRANSFERASE TEST R1	R1 4×500ml
BC-0007G	MeDiPro ALANINE AMINOTRANSFERASE TEST R2	R2 4×200ml

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan