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台塑澱粉酶試劑（AMY）- colorimetric test-CNPG3 method

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MeDiPro α-AMYLASE TEST (AMY) - colorimetric test-CNPG3 method

INTENDED USE


For the quantitative determination of α-amylase activity in serum, plasma or urine.

CLINICAL SIGNIFICANCE

α-Amylase catalyzes the hydrolysis of oligosaccharides, polysaccharides, starches and glycogen with formation of maltose and dextrin. The human serum α-amylase activity is therefore clinically related to the diseases of the pancreas and for evaluation of pancreatic function. α-amylase activity increased considerably in acute pancreatitis and obstruction of the pancreatic ducts. α-amylase activity also increase in acute abdominal pain associated with peptic ulcer, empyema of the gallbladder and intestinal obstruction.

PRINCIPLE

Substrate Gal-G₂-α-CNP is hydrolyzed by α-amylase to Gal-G₂ and CNP stoichiometrically. The rate of CNP formation due to substrate hydrolysis by α-amylase is proportionally correlated with α-amylase activity which is measured by following the rate of absorbance increase at 410nm. The reaction is shown as follows:

REAGENT

- Package: please see the reagent box label shown.
- Components:

Component	Conc.
Gal-G ₂ -α-CNP	0.3 mmol/L
Buffer	50 mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at 2 – 8 °C., reagent should be protected from light and contamination should be avoided. **Do not freeze the reagent!**

SPECIMEN COLLECTION AND PREPARATION

Serum, heparin-treated plasma, or urine samples may be used. Anticoagulant such as citrate, oxalate or EDTA must be avoided. α-Amylase activity in serum or urine sample, without bacterial contamination, are stable for 7 days at room temperature and several months at 4°C.

PROCEDURES

- Main wavelength : 405 nm
Sub. wavelength : 660nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains single reagent, ready to use.

	Blank	Control	Specimen
Specimen ml	---	---	0.02
Control ml	---	0.02	---
ddH ₂ O ml	0.02	---	---
Reagent ml	1.0	1.0	1.0

Mix, incubate at 37°C for 5 min, and read the initial

absorbance A₁ against reagent blank, then read end absorbance A₂ in every 1 min. A = A₄₀₅ - A₆₆₀

CALCULATION

$$\text{AMY (U/L)} = \frac{(A_2 - A_1) / \text{min} \times V_t \times 1000}{L_p \times \epsilon \times V_s}$$

$$= (A_2 - A_1) / \text{min} \times 3806$$

V_t: Reaction total volume 1.02 ml, V_s: sample volume 0.02 ml
ε: CNP molar absorptivity 13.4,
1000: transfer U/ml to U/L, L_p: Optical path length (cm)

REFERENCE RANGE

Serum < 220 U/L
Fresh urine < 1000 U/L
Collect urine < 900 U / 24 h

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition.
- The test is developed to determine α-amylase concentrations up to 2000U/L. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- Anticoagulant such as citrate, oxalate or EDTA must be avoided because of binding of calcium ion, which is essential for α-amylase activity. Specimen lipemia might be affect test result. Saliva contains very high α-amylase activity; avoid bringing any saliva into the substrate.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance

MeDiPro α-AMYLASE TEST (AMY) - colorimetric test-CNPG3 method

with National or local regulations related to the safety precautions of such materials.


- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)
- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.
- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[U/L]	SD [U/L]	CV[%]
Sample1	57.90	1.29	2.22
Sample2	384.00	2.91	0.76
Sample3	380.40	3.06	0.81

- Linearity

This kit has a good linearity up to 3000U/L.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 500mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 0.2%

- Stability

Expire day	1 year
Open vial stability	14 day

REFERENCE

- Winn-Deen E.S., David H., Sigler G., Chavez R. Development of a direct assay for α-amylase. 1988. Clin. Chem. 34:2005-2008.
- Foo A.Y., Bias R. Amylase measurement with 2-chloro-4-nitrophenyl maltotriose as substrate. 1998. Clin. Chim. Acta. 272:137-147.

PARAMETER SETUP

Hitachi 7170/917 Applications

TEST	[AMY]
ASSAY CODE	[RateA]:[3]-[10]
SAMPLE VOLUME	[4]
R1 VOLUME	[200]
R2 VOLUME	[0]
WAVELENGTH(nm)	[660][405]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[AMY]
ASSAY CODE	[RateA]:[5]-[15]
SAMPLE VOLUME	[6]
R1 VOLUME	[300]
R2 VOLUME	[0]
WAVELENGTH(nm)	[660][405]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0009M	MeDiPro α-AMYLASE TEST	R1 6×20ml
BC-0009A	MeDiPro α-AMYLASE TEST	R1 4×60ml
BC-0009B	MeDiPro α-AMYLASE TEST	R1 4×100ml
BC-0009C	MeDiPro α-AMYLASE TEST	R1 2×250ml
BC-0009D	MeDiPro α-AMYLASE TEST	R1 2×500ml

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan