### 台塑血中尿素氮試劑 (BUN) - Enzymatic UV test-urease-GLDH

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## MeDiPro BLOOD UREA NITROGEN TEST (BUN) - Bromcresol green method

### INTENDED USE


For the quantitative determination of urea nitrogen in serum or plasma

### CLINICAL SIGNIFICANCE

Urea is a nitrogen-generated end-product of proteolysis. Urea is produced in the liver. It constitutes the majority of the NPN fraction of the blood, and is normally excreted by the kidney into urine. Blood Urea Nitrogen (BUN) levels are therefore related to the protein metabolism (Intake and catabolism) and kidney function.

### PRINCIPLE

Urea is hydrolyzed by urease to produce ammonia. The ammonia is then coupled with  $\alpha$ -ketoglutarate and NADH to produce glutamate and NAD<sup>+</sup>. The rate of absorbance decrease is directly proportional to the amount of urea present in the sample. The reaction equation is as following:

### REAGENT

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R <sub>1</sub> :	Tris buffer	100 mmol/L
	$\alpha$ -ketoglutarate	5 mmol/L
	ADP	0.6 mmol/L
	GLDH	>1500 U/L
	Urease	>20 U/L
R <sub>2</sub> :	Tris buffer	100 mmol/L
	NADH	0.25 mmol/L

### STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at 2 – 8 °C., reagent should be protected from light and contamination should be avoided.

**Do not freeze the reagent!**

### SPECIMEN COLLECTION AND PREPARATION

The test can be performed by using serum, plasma. For serum, blood is drawn into a tube which does not contain anticoagulant and allowed to form clotting. The serum is then separated from the clot. Samples are generally stable for 1 day at room temperature but can be stored for several days refrigerated or several months frozen without any appreciable loss of urea.

For plasma, add whole blood directly into a tube containing anticoagulant. Most common anticoagulants may be used except EDTA.

### PROCEDURES

- Main wavelength : 340nm  
Sub. wavelength : 405nm  
Reaction Temperature : 37°C  
Optical path length : 1.0 cm
- This kit contains two reagents, ready to use.

	Control	Specimen
R <sub>1</sub> ( ml )	0.8	0.8
Control ml )	0.01	----
Specimen ( ml )	----	0.01
Mix, 37°C incubate 5 min		
R <sub>2</sub> ( ml )	0.2	0.2

Mix, incubate at 37°C for 1 min, and read the initial absorbance A<sub>1</sub> against reagent blank, then read end absorbance A<sub>2</sub> in every 30 sec. A = A<sub>340</sub>-A<sub>405</sub>.

### CALCULATION

With standard or calibrator

$$\text{BUN (mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

### REFERENCE RANGE

5~25 mg/dL (1.8~8.9 mmol/L)

### WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition.
- Each instrument, material and distilled water used in the test procedures shall be confirmed without ammonium contamination.
- The test is developed to determine blood urea nitrogen concentrations up to 150mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- Do not use if reagent blank less than 1.0.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in

## MeDiPro BLOOD UREA NITROGEN TEST (BUN) - Bromcresol green method

Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.


- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)
- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.
- Validity please see the reagent box label shown.

### REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	15.50	0.26	1.71
Sample2	48.29	0.29	0.61
Sample3	49.43	0.49	1.00

- Linearity

This kit has a good linearity up to 150mg/dL.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 500mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 1.8%

- Stability

Expire day	1 year
Open vial stability	30 day

### REFERENCE

Tiffany, TO, Jansen JM, Burtis CA, et.al. Enzymatic kinetic rate and endpoint analysis of substrate by use of GEMSAEC fast analyzer. Clin. Chem. 1972. 18:829.

### PARAMETER SETUP

#### Hitachi 7170/917 Applications

TEST	[BUN]
ASSAY CODE	[Rate A]: [19]-[22]
SAMPLE VOLUME	[2]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH(nm)	[405][340]
CALIB. METHOD	[Linear]

#### Hitachi 7150/717 Applications

TEST	[BUN]
ASSAY CODE	[Rate A]: [30]-[35]
SAMPLE VOLUME	[3]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH(nm)	[405][340]
CALIB. METHOD	[Linear]

### ORDERING INFORMATION

Cat. No.	Product	Package
BC-0012M	MeDiPro BLOOD UREA NITROGEN TEST	R1 6×20ml R2 3×10ml
BC-0012A	MeDiPro BLOOD UREA NITROGEN TEST	R1 4×60ml R2 2×30ml
BC-0012B	MeDiPro BLOOD UREA NITROGEN TEST	R1 4×100ml R2 2×50ml
BC-0012C	MeDiPro BLOOD UREA NITROGEN TEST R1	R1 2×300ml
BC-0012D	MeDiPro BLOOD UREA NITROGEN TEST R1	R1 2×500ml
BC-0012G	MeDiPro BLOOD UREA NITROGEN TEST R2	R2 2×200ml

### FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

### FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>  
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381  
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan