台塑直接膽紅素試劑（DBIL）- Diazo-DPD method

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MeDiPro DIRECT BILIRUBIN TEST (DBIL) - Diazo-DPD method


INTENDED USE

For the quantitative determination of direct bilirubin in serum or plasma.

CLINICAL SIGNIFICANCE

An increase in the formation or retention of bilirubin in the body results in the increased level of serum bilirubin and jaundice. This hyperbilirubinemia is classified as either pre-hepatic, hepatic or post-hepatic depending on the major cause of condition. Therefore, the determination of the total bilirubin and its conjugated (direct) bilirubin is important for the diagnosis of hyperbilirubinemia.

PRINCIPLE

REAGENT

- Package: please see the reagent box label shown.
- Components:

	Component	Conc.
R ₁	Sulfamic acid	70 mmol/L
	HCl	70 mmol/L
R ₂	3,5-DCPD	5mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at **2 – 8 °C.**, reagent should be protected from light and contamination should be avoided.
Do not freeze the reagent!

SPECIMEN COLLECTION AND PREPARATION

Fresh serum and heparin-treated plasma are the choices. The sample should be collected without hemolysis, since hemoglobin will inhibit the diazo reaction. Avoid direct light exposure to the specimen since bilirubin values may decrease as much as 50% in one hour. Serum specimen may be kept in dark in 2-8°C for up to one week, and in freezer for 3 months without appreciable change in the bilirubin levels.

PROCEDURES

- Main wavelength : 546 nm
Sub. wavelength : 660nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains two reagents, ready to use.

	Blank	Control	Specimen
R ₁ (ml)	0.8	0.8	0.8
ddH ₂ O (ml)	0.04	---	---
control (ml)	---	0.04	---
Specimen(ml)	---	---	0.04
Max, 37°C incynate 5min			
R ₂ (ml)	0.2	0.2	0.2

Mix, incubate at 37°C for 5 min, and read the absorbance against reagent blank. $A = A_{546} - A_{660}$.

CALCULATION

With standard or calibrator

$$\text{DBIL (mg/dL)} = \frac{A_{\text{sample}}}{A_{\text{std./cali.}}} \times \text{conc. Std./cali. (mg/dL)}$$

REFERENCE RANGE

0.2-0.7 mg/dL

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Hemolysis and lipemia specimen should be avoided.
- The test is developed to determine direct bilirubin concentrations up to 6mg/dL. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan

MeDiPro DIRECT BILIRUBIN TEST (DBIL) - Diazo-DPD method


in accordance with National or local regulations related to the safety precautions of such materials.

- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on potentially hazardous components. (MSDS could be obtained from local dealer.)
- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.
- Validity please see the reagent box label shown.

REAGENT CHARACTERS

1. Precision (Within run)			
N=15	Mean[mg/dL]	SD [mg/dL]	CV[%]
Sample1	0.54	0.01	2.27
Sample2	2.48	0.02	0.71
Sample3	2.36	0.03	1.45

- Linearity

This kit has a good linearity up to 12mg/dL.

- Interference

Interference	Influence effect
Hemoglobulin	Not suitable when hemolysis occur
Ascorbic acid	No interference was observed by ascorbic acid up to 35mg/dL
Intrafat	Not suitable when lipemia occur

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

- Doumas BT, Wu TW. The measurement of bilirubin fractions in serum. Crit Rev Clin Lab Sci 1991; 28(5-6): 415-445.
- Billy W.P., Basll T.D., etc. A candidate reference method for determination of bilirubin in serum. Test for transferability. Clin. Chem. 1983, 29: 297-301

PARAMETER SETUP

Hitachi 7170/917 Applications

TEST	[DBIL]
ASSAY CODE	[2POINT]: [16]-[34]
SAMPLE VOLUME	[8]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH(nm)	[700][570]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[DBIL]
ASSAY CODE	[2POINT]: [24]-[50]
SAMPLE VOLUME	[12]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH(nm)	[700][570]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0011M	MeDiPro DIRECT BILIRUBIN	R1 6×20ml
	TEST	R2 3×10ml
BC-0011A	MeDiPro DIRECT BILIRUBIN	R1 4×60ml
	TEST	R2 2×30ml
BC-0011B	MeDiPro DIRECT BILIRUBIN	R1 4×100ml
	TEST	R2 2×50ml
BC-0011C	MeDiPro DIRECT BILIRUBIN	R1 2×300ml
	TEST R1	
BC-0011D	MeDiPro DIRECT BILIRUBIN	R1 2×500ml
	TEST R1	
BC-0011G	MeDiPro DIRECT BILIRUBIN	R2 2×200ml
	TEST R2	

FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan