



MeDiPro γ-GLUTAMYL TRANSFERASE TEST (GGT) - colorimetric test


INTENDED USE

For the quantitative determination of γ-glutamyl transferase activity in serum or plasma.

CLINICAL SIGNIFICANCE

γ-glutamyl transferase (GGT) is a membrane-localized enzymes that catalyzes the transfer of a γ-glutamyl group from a γ-glutamylpeptide to another peptide or amino acid. Kidney, pancreas and liver are rich in γ-GT. Serum γ-GT is generally elevated as a result of liver disease. Cholestasis caused by alcohol or drug ingestion, mechanical or viral cholestasis, liver metastases will increase the activity of γ-GT. In bone disorders which alkaline phosphatase is elevated but γ-GT is normal; and in skeletal muscle disorder which the AST is elevated but γ-GT is normal.

PRINCIPLE



REAGENT

- Package: please see the reagent box label shown.
- Components:

Component	Conc.
R ₁ : Tris buffer	110 mmol/L
glycylglycine	100 mmol/L
R ₂ : L-γ-glutamyl-p-nitroanilide	25 mmol/L

STORE TEMPERATURE

The standard is stable up to the end of the indicated expiration date. If stored at 2 – 8 °C., reagent should be protected from light and contamination should be avoided. **Do not freeze the reagent!**

SPECIMEN COLLECTION AND PREPARATION

Serum or lithium heparized plasma is the choice. EDTA and citrate do not interfere with analysis. Serum GGT values are stable at room temperature or 4°C for at least 7 days and are stable for at least 2 months when frozen.

PROCEDURES

- Main wavelength : 405nm
Sub. wavelength : 660nm
Reaction Temperature : 37°C
Optical path length : 1.0 cm
- This kit contains two reagents, ready to use.

	Blank	Control	Specimen
R ₁ (ml)	0.8	0.8	0.8
ddH ₂ O (ml)	0.05	---	---
control (ml)	---	0.05	---
Specimen(ml)	---	---	0.05
R ₂ (ml)	0.2	0.2	0.2

Mix, incubate at 37°C for 5 min, and read the absorbance against reagent blank. $A = A_{405} - A_{660}$

CALCULATION

$$\text{GGT (U/L)} = \frac{(A_2 - A_1) / \text{min} \times V_t \times 1000}{L_p \times \epsilon \times V_s} = (A_2 - A_1) / \text{min} \times 2210$$

V_t: Reaction total volume 1.05 ml, V_s: sample volume 0.05 ml
ε: p-nitroaniline molar absorptivity 9.5,
1000: transfer U/ml to U/L, L_p: Optical path length (cm)

REFERENCE RANGE

Sex	Conc. (U/L)
Male	11-50 U/L
Female	7-32 U/L

WARNINGS AND PRECAUTIONS

- This kit offers an optional calibrator, which is sold individually. Bio-Rad Lyphochek control is recommended to use as serum control.
- Each laboratory has to perform the quality control test to assure the results being reliable before running the specimen tests.
- This kit is for professionals and *in vitro* diagnostic use only.
- To ensure the accuracy of result, the absorbance should be measured within 30 minutes after sample addition. Hemolysis and lipemia might interfere with the test result.
- The test is developed to determine γ-glutamyl transferase concentrations up to 800U/L. When values exceed this range, samples should be diluted with normal saline and calculate the results by multiplying the dilution factor.
- Do not use if reagent turbid or blank OD higher than 1.0.
- The above-mentioned procedures are suitable either for the general semi-automatic, full-automatic biochemical analysis instrument or manual operation.
- Since all specimens are potentially infectious, they should be handled with appropriate precautions and practices in accordance with Biosafety level 2 as recommended by USA NIH manual Biosafety in Microbiological and Biomedical Laboratories, and in accordance with National or local regulations related to the safety precautions of such materials.
- Waste management please refers to the local legal requirements.
- Please refer to the manufacturer's safety data sheet and the product labeling for information on



MeDiPro γ-GLUTAMYL TRANSFERASE TEST (GGT) - colorimetric test

potentially hazardous components. (MSDS could be obtained from local dealer.)


- According to the technical suggestion, the volume of reagent and specimen could be adjusted in a ratio for full-automatic biochemical analysis instrument use. It won't affect the absorbance and the result.
- Validity please see the reagent box label shown.

REAGENT CHARACTERS

- Precision (Within run)

N=15	Mean[U/L]	SD [U/L]	CV[%]
Sample1	35	0.85	2.46
Sample2	191	0.82	0.43
Sample3	197	0.84	0.43

- Linearity



This kit has a good linearity up to 1200U/L.

- Interference

Interference	Influence effect
Hemoglobin	No interference was observed by hemoglobin up to 150mg/dL
Ascorbic acid	No interference was observed by ascorbic acid up to 50mg/dL
Bilirubin (free form)	No interference was observed by bilirubin up to 40mg/dL
Bilirubin (conjugate form)	No interference was observed by bilirubin up to 40mg/dL
Intrafat	No interference was observed by intrafat up to 0.4%

- Stability

Expire day	1 year
Open vial stability	30 day

REFERENCE

- IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C, Part 6. Clin. Chem. Lab. Med. 2002, 40: 631.
- New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clinica Chimica. Acta. 2003, 327: 69.

PARAMETER SETUP

Hitachi 7170/917 Applications

TEST	[GGT]
ASSAY CODE	[Rate A]: [19]-[34]
SAMPLE VOLUME	[10]
R1 VOLUME	[160]
R2 VOLUME	[40]
WAVELENGTH (nm)	[660][405]
CALIB. METHOD	[Linear]

Hitachi 7150/717 Applications

TEST	[GGT]
ASSAY CODE	[Rate A]: [30]-[50]
SAMPLE VOLUME	[15]
R1 VOLUME	[240]
R2 VOLUME	[60]
WAVELENGTH(nm)	[660][405]
CALIB. METHOD	[Linear]

ORDERING INFORMATION

Cat. No.	Product	Package
BC-0018M	MeDiPro γ-GLUTAMYL TRANSFERASE TEST	R1 6×20ml R2 3×10ml
BC-0018A	MeDiPro γ-GLUTAMYL TRANSFERASE TEST	R1 4×60ml R2 2×30ml
BC-0018B	MeDiPro γ-GLUTAMYL TRANSFERASE TEST	R1 4×100ml R2 2×50ml
BC-0018C	MeDiPro γ-GLUTAMYL TRANSFERASE TEST	R1 2×300ml
BC-0018D	MeDiPro γ-GLUTAMYL TRANSFERASE TEST	R1 2×500ml
BC-0018G	MeDiPro γ-GLUTAMYL TRANSFERASE TEST	R2 2×200ml



FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan



FORMOSA BIOMEDICAL TECHNOLOGY CORP.

F-5F, No. 201, Tunghua N. Rd, Taipei, 105, Taiwan Website: <http://www.fbc.com.tw/>
TEL: +886-2-2712-2211 #7822 FAX: +886-2-2717-8381
Factory: No. 3, Longchuan Rd, Longtang Village, Jiaosi, Yilan County, 262, Taiwan